

Variability of fungi detection in the respiratory tract of racehorses according to the sampling site

Louise C. LEMONNIER¹, A. COUROUCE^{1,2}, M. CESSANS¹, L. PETIT³, J.M. CARDWELL⁴, M-P. TOQUET³, E.A. RICHARD^{2,3}

¹Department of Equine Internal Medicine, University Animal Hospital, CISCO-Oniris, Route de Gachet, 44300 Nantes, France

²BIOTARGEN, Normandie Univ, UNICAEN, 14000 Caen, France

³Labeo-Frank Duncombe Laboratory, 1 route de Rosel, 14053 Caen, France

⁴Department of Pathobiology and Population Sciences, Royal Veterinary College, Hatfiled, Hertfordshire, UK





INTRODUCTION

- The clinical relevance of colonization of equine airways by fungi is unclear.
- Diagnosis of fungi in the respiratory tract is fastidious and there is no method of reference.

OBJECTIVES

- To characterize fungi from upper and lower airways of racehorses
- To determine the most reliable methodology for fungal detection

HYPOTHESIS

• Fungal cytology of a pooled BAL (BAL_P) is representative of corresponding left (BAL_L) and right (BAL_R) samples.

METHODS

- * 30 Standardbred racehorses in active training, aged 4.6 \pm 1.6 years, from 3 different training stables
- TW, followed by BAL from each lung separately, collected through a guarded catheter and gastroscope
- Fungi detection, identification and characterization by cytology and culture
 - Positive culture: ≥ 1 mold identified
 - Positive cytology: ≥ 1 fungal element (hyphae, spore) observed

RESULTS

- Fungal culture (Table 1)
 - <u>TW</u>: 28/30 (93%) positive
 - BAL_E: 4/30 (13%) positive → All positive in TW, 3/4 with same fungi
 - 9 genus (alone or concomitant): Aspergillus (5 genus)
 13/28 (46%), Penicillium 12/28 (43%)

Table 2: Cytology of pooled and individual BALF, and of TW

	Cytology positive	Spores	Hyphae	Phagocyted
TW	22/30	18/22	16/22	14/22
BAL _P	7/30 3 bilateral 4 unilateral	7/7	2/7	6/7
BAL _R	8/30	5/8	6/8	5/8
BAL	7/30	7/7	0/7	6/7

- Comparison of sample sites (Table 3)
 - Poor agreement between BAL_I and BAL_R
 - Poor agreement between culture and cytology
 - BAL_P vs separate lungs:
 - Not equivalent (p = 0.025 by McNemar Test)
 - Sensitivity = 0,58; Specificity = 1,00

<u>Table 1</u>: Identification of fungi in horses with a positive BAL_P culture

	Case number		9		17		18		21
	Fungi genera	TW	BAL _P	TW	BAL₽	TW	BAL₽	TW	BAL _P
	A fumigatus							Х	
Ш	A nidulans	х		х	х	Х			
ľ	A versicolor							Х	
П	Penicillium	х		х		Х	Х	Х	х
I	Aureobasidium		Х						
	Chrysosporum	Х							

Fungal cytology (Table 2)

- <u>TW</u>: 22/30 (73%) positive; 12/22 (55%) with both spores and hyphae, frequently phagocyted
- BAL: Large disparity between lungs

 $\underline{\textit{Table 3}} : \textit{Comparison of BALF sample sites, and of BALF methods of analysis}$

	BAL _L only	BAL _R only	BAL _{L+R}	Agreement (κ)	
Cytology	4/7	5/8	3	0,015	
	Culture only	Cytology only	Culture + cytology	Agreement (κ)	

DISCUSSION

- Despite a similar prevalence of fungi detection in left and right BAL, a large disparity was found between lungs on cytology.
- Detecting fungi by cytology on either BAL_P or BAL_L + BAL_R was not equivalent

CONCLUSION

Pooled BAL is a highly specific while fairly sensitive alternative to combined cytology of individual left/right samples for mold/fungi detection.







