

#### **BLEPal**

Mohammed Farouk Akli, Nedir Nour El Amine Boukerras, Nesrine Derradji, Dyhia Laga, Mohammad Imran Syed

#### ▶ To cite this version:

Mohammed Farouk Akli, Nedir Nour El Amine Boukerras, Nesrine Derradji, Dyhia Laga, Mohammad Imran Syed. BLEPal: Bluetooth Trace Synchronization and Merging Python Tool. [Technical Report] Sorbonne Université UPMC. 2022. hal-03765103v1

## HAL Id: hal-03765103 https://hal.science/hal-03765103v1

Submitted on 30 Aug 2022 (v1), last revised 19 Sep 2022 (v2)

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

# **BluePal**

# Bluetooth Trace Synchronization and Merging Python Tool

Mohammed Farouk AKLI, Nedir Nour El Amine BOUKERRAS, Nesrine DERRADJI, Dyhia LAGA

(All authors contributed equally to this work)

Supervisor: Mohammad Imran SYED



# **Table of Contents**

1. Introduction		oduction	3
	1.1.	Synchronization Process	3
	1.2.	Merging Process	5
2.	How to use the tool		6
	2.1.	Libraries required	6
	2.2.	The configuration file file.cfg	6
3.	Link	x for the tool	7
Re	References		



1. Introduction 3

#### 1. Introduction

The application we developed is written in python and shell script, it encompasses every aspect of the process: from the capture procedure until the final result. Every capture result will be written on a .csv file, the user would not need to do anything manually. Our application is based on an automation philosophy (everything is happening in an automatic way transparent to the end user) and a user friendly dialog inbox which will accompany the user for the whole process.

We decided to name it BluePal, this name is inspired by the tool WiPal [1] that was made for Wi-Fi, as BluePal is a tool made for Bluetooth and pal means a friend so our application is metaphorically a Bluetooth's friend as it offers a trace completeness for the protocol so we can higher the performances and reduce data loss when applying the synchronizing and merging process.

To use BluePal, you will only need the python script of the application and a configuration file named file.cfg. The python script of BluePal and the file.cfg need to be in the same directory.

The application uses either the nrf52840 dongle or the Ubertooth one dongle to capture the Bluetooth traffic.

#### 1.1. Synchronization Process

Our synchronization method is based on synchronizing an asynchronous environment.

Our environment is a set of Raspberry Pi 4 devices and Bluetooth sniffers (either **Nordic** or **Ubertooth**), a Bluetooth headset and a PC. we establish our topology as following (caricatured by a graph):




Figure 1.1: Capture's graph



1. Introduction 4

A timestamp is associated with every **Raspberry** + **sniffer** set, the timestamp of the sniffer will be equal to the timestamp of the Raspberry, we would not assume that we are in a synchronous environment (every sniffer can have its own timestamp).

Let us associate a timestamp to every sniffer of the four sniffers that we have at the exact moment before launching the capture:

Timestamp of sniffer 1 :  $t_1$ 

Timestamp of sniffer 2:  $t_2$ 

Timestamp of sniffer 3:  $t_3$ 

Timestamp of sniffer 4:  $t_4$ 

After launching the capture and the Bluetooth traffic, a specific frame has its own timestamp (the instant the frame was captured by the sniffer). This timestamp will be, obviously, bigger than the initial timestamp before launching the capture. Let us associate a timestamp to a frame of every sniffer of the four sniffers that we have (for simplification we considered here a single frame for every sniffer but in reality we have many frames within a sniffer's trace):

Timestamp of the frame of the sniffer 1:  $t_1 + x_1$ 

Timestamp of the frame of the sniffer 2:  $t_2 + x_2$ 

Timestamp of the frame of the sniffer 3:  $t_3 + x_3$ 

Timestamp of the frame of the sniffer 4:  $t_4 + x_4$ 

Now and by doing a subtraction between the two set of equations we will have a new timestamp, :

New frame's timestamp of the sniffer 1:  $t_1 + x_1 - t_1 = x_1$ 

New frame's timestamp of the sniffer 2:  $t_2 + x_2 - t_2 = x_2$ 

New frame's timestamp of the sniffer 3:  $t_3 + x_3 - t_3 = x_3$ 

New frame's timestamp of the sniffer 4:  $t_4 + x_4 - t_4 = x_4$ 

These new timestamps  $x_i$  are absolute timestamps rather than the relative timestamps we had before the subtraction operation, they do not depend on the timestamp of the machine (the timestamp of the sniffer).

We can now establish an order to these absolute timestamps, for example we can have this following order:

$$x_1 < x_2 < x_3 < x_4$$

So the frame with timestamp  $x_1$  is the first captured frame of the experiment, the frame with the timestamp  $x_2$  is the second frame captured, the frame with the timestamp  $x_3$  is the third frame captured and the frame with the timestamp  $x_4$  is the last frame captured.



1. Introduction 5

#### 1.2. Merging Process

The merging process will lean on the extraction of unique frames. When the capture scenario uses **Nordic** sniffers, we will remove the duplicates according to the couple [CRC,Channel Index]. When the capture scenario uses **Ubertooth** sniffers, we will remove the duplicates according to the value of the CRC.

Here is an example of how the merging process works (example given for 4 traces):




Figure 1.2: Merging process

The final result will be the trace containing all four traces synchronized and merged (trace 1, 2, 3 and 4).



2. How to use the tool 6

#### 2. How to use the tool

On a console, type python3 bluepal.py then a dialog inbox will guide you through the whole process<sup>1</sup>.

#### 2.1. Libraries required

You need to have the following libraries installed:

- os
- pandas
- subprocess
- re
- signal
- time
- matplotlib.ticker
- tkinter
- sys
- csv

#### 2.2. The configuration file file.cfg

the content of file.cfg:

<sup>1.</sup> It is necessary to use Python3.



3. Link for the too

## 3. Link for the tool

The tool can be found on the following link: <a href="https://gitlab.lip6.fr/syed/bluepal">https://gitlab.lip6.fr/syed/bluepal</a>



References 8

#### References

[1] Thomas Claveirole and Marcelo Dias de Amorim. 2010. WiPal: efficient offline merging of IEEE 802.11 traces. SIGMOBILE Mob. Comput. Commun. Rev. 13, 4 (March 2010), 39–46. DOI: <a href="https://doi.org/10.1145/1740437.1740443">https://doi.org/10.1145/1740437.1740443</a>

