

DTA characterization of TBPB and TBPB-CO 2 hydrates in porous media

Pascal Clain, Anthony Delahaye, Laurence Fournaison

▶ To cite this version:

Pascal Clain, Anthony Delahaye, Laurence Fournaison. DTA characterization of TBPB and TBPB-CO 2 hydrates in porous media. CTEC, Jun 2012, Lyon, France. hal-03781506

HAL Id: hal-03781506 https://hal.science/hal-03781506

Submitted on 26 Sep 2022

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



DSC characterization for TBPB hydrate in porous media



Pascal CLAIN*, Anthony DELAHAYE, Laurence FOURNAISON LGP2ES (EA21), Irstea GPAN, 1 rue Pierre-Gilles de Gennes, CS 10030 92761 Antony cedex



Sensors of T. P

Thermal studies

Objectives of this study

Coupling between a storage device and a multi temperature cold system

Modeling/Conception of an unique hydrate storage at various temperature levels

Interest of hydrate slurry for refrigeration

- Crystals with high melting latent heat: up to ice
- $T_{eq} > 0^{\circ}$ suitable for air conditioning and refrigeration
- No mechanical generation

Previous studies in Irstea

- CO₂+additives
- Rheology hydrate slurries in flow

Equilibrium conditions of TBPB hydrate in porous media

The method of extrapolation temperature is used for determine the dissociation temperature at atmospheric pressure on DSC curves



Dissociation temperature in porous media is lower than bulk media. At stoichiometric concentration, pore sizes have not effect



DSC device

Determination of solid mass hydrate fraction

Experimental set-up: 1DSC and 1 DTA device

Solid hydrate fraction calculated from a mass balance on CO₂

Hydrate mass fraction is determined by a calorimetric method, ba

Hydrate mass fraction is determined by a calorimetric method, based on the composition of TBPB hydrate.

DTA device

$$x_{h} = \frac{\Delta H_{total} - \Delta H_{ice}}{\Delta H_{ice} - \Delta H_{hydrate}}$$

Porous media are a high specific surface. The formation of hydrates in porous media have to improve. The results show that the amount formed in porous media is equivalent to bulk media although the formation is more difficult.



Porous media is a hydrate inhibitor . Due to the capillary effects it is more difficult to form hydrates in pore. Dissociation temperature shift down in porous media. At a given concentration of salt and given pressure, melting temperature is lowest when the pore radius is minimal. That could be an interesting way for the multi temperature storage Amount of hydrate formed in porous media is equivalent to bulk phase. For each concentration, maximum hydrate formed is obtained at the largest pore size PERSPECTIVES - Currently, a thermal analysis is in progress to complete the characterization of TBPB in presence of CO₂ on DTA device

Acknowledgement : This work is supported by ADEME