



# High temperature oxidation behaviour of alloys containing very high densities of carbides

Patrice Berthod, Pierrick Lemoine, Lionel Aranda

## ► To cite this version:

Patrice Berthod, Pierrick Lemoine, Lionel Aranda. High temperature oxidation behaviour of alloys containing very high densities of carbides. Eurocorr 2007, Sep 2007, Freiburg im Breisgau, Germany.  
. hal-03906066

HAL Id: hal-03906066

<https://hal.science/hal-03906066>

Submitted on 19 Dec 2022

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

# High temperature oxidation behaviour of alloys containing very high densities of carbides


Patrice BERTHOD, Pierrick LEMOINE, Lionel ARANDA

Laboratoire de Chimie du Solide Minéral, Université Henri Poincaré, Nancy-Université  
BP 239, Bd des Aiguillettes, 54506 Vandoeuvre lès Nancy, France


Some alloys containing both very high fractions of carbides and a high amount of chromium can potentially display a low thermal expansion and a good resistance against high temperature oxidation, simultaneously. The aim of this work is to explore the high temperature oxidation behaviour of such alloys, by studying three alloys based on nickel and one based on cobalt, all of them containing more than one weight percent of carbon to obtain very high carbides fractions and 30%wt. of chromium added to achieve a chromia-forming behaviour.

Four alloys were prepared by induction melting in a cold crucible under inert atmosphere (300mbar of pure Argon). Three are based on Ni-30%Cr with three different high carbon contents (1.2%C, 1.6%C and 2.0%C, weight percent) and a fourth one is based on Co-30%Cr with a high fraction of tantalum carbides achieved by 0.8%C and 12%Ta (for equal atomic contents between carbon and tantalum). Thermogravimetry tests were performed at 1000°C, 1100°C and 1200°C for 50 hours in dry air.


Alloy	Ni	Co	C	Cr	Ta
Ni12	Bal.	/	1.2	30	/
Ni16	Bal.	/	1.6	30	/
Ni20	Bal.	/	2.0	30	/
CoTa08	/	Bal.	0.8	30	12



Mass gain curves obtained for all the alloys when oxidized at 1200°C (top), 1100°C (middle) and 1000°C (bottom)



{dm/dt = Kp - Kv . m}-type treatment of mass gain results and validation of the obtained Kp and Kv by comparison between experimental curve and the mathematical curve plotted with these values (here: CoTa08 at 1200°C)



Values of Kp and Kv issued from the {dm/dt = Kp - Kv . m} treatment of the thermogravimetry results for all alloys at the three temperatures

Kp ( $\times 10^{-12}$ $\text{g}^2 \cdot \text{cm}^{-4} \cdot \text{s}^{-1}$ )	temperature	Ni12	Ni16	Ni20	CoTa08
		Ni30Cr 1.2%C	Ni30Cr 1.6%C	Ni30Cr 2.0%C	Co30Cr 0.8%C 12%Ta
		1200°C	103	65	344
Kv ( $\times 10^{-10}$ $\text{g} \cdot \text{cm}^{-2} \cdot \text{s}^{-1}$ )	temperature	Ni12	Ni16	Ni20	CoTa08
		1200°C	102	117	/
		1100°C	59	85	94
		1000°C	27	24	7.6

**Bold characters:** excellent fit between the experimental mass gain curve and the mathematical one  
**Normal characters:** good fit between the experimental mass gain curve and the mathematical one

The studied alloys are almost all sufficiently resistant against oxidation for this exposure duration. The parabolic constant Kp and the constant of chromia volatilization Kv for the nickel alloys do not follow a systematic law of variation versus the carbon content. However the depth of the carbide-free zone clearly decreases when the carbon content increases. It is possible that the impoverishment of the matrix for very high fractions of chromium carbides (which then contain a great part of chromium), can lead to local losses of the chromia-forming behaviour.

In the case of the cobalt alloy the presence of tantalum obviously leads to fast mass gain kinetics but delays catastrophic oxidation.