

Differentiation of the ant genus Tapinoma from the Mediterranean Basin by species-specific cuticular hydrocarbon profiles.

Laurence Berville, Abraham Hefetz, Alain Lenoir, Marielle Renucci, O Blight, Provost E.

▶ To cite this version:

Laurence Berville, Abraham Hefetz, Alain Lenoir, Marielle Renucci, O Blight, et al.. Differentiation of the ant genus Tapinoma from the Mediterranean Basin by species-specific cuticular hydrocarbon profiles.. THE 12TH INTERNATIONAL COURSE IN CHEMICAL ECOLOGY, Nov 2012, Iena, Germany. Unpublished, 2012, 10.13140/RG.2.2.27178.39360. hal-04212281

HAL Id: hal-04212281

https://hal.science/hal-04212281

Submitted on 20 Sep 2023

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.













Differentiation of the ant genus Tapinoma from the Mediterranean Basin by species-specific cuticular hydrocarbon profiles


Correct species identification is a precondition for biological study; yet despite a long history of morphological investigations, the systematic position of many ant species remains unclear. Like any phenotypic character, cuticular hydrocarbons (CHCs) are reliable indicators of species identity. Here, we compared and identified cuticular hydrocarbon profiles of workers of several species of *Tapinoma*.


50 Colony fragments were collected in **Algeria**, **France** (Corsica and continental France), **Morocco**, **Israel** and **Spain**. The nesting sites sampled ranged over **1900** km from the northern to the southern part and over **4800** km from the eastern to the western part of the study area.

Berville L^{1*}., Hefetz A²., Lenoir A³., Renucci M¹., Blight O¹., & E. Provost¹

¹ Aix-Marseille Université, CNRS, IMBE, 7263 UMR, Campus Aix Technopôle Arbois-Méditerranée, Pavillon Villemin, F-13545 Aix-en-Provence Cedex 04, France. ² Department of Zoology, George S. Wise Faculty of Life Sciences, Tel Aviv University, Ramat Aviv, 69978, Israel. ³ IRBI, UMR CNRS 7261, Université François Rabelais, Faculté des Sciences et Techniques, Parc Grandmont, 37200 Tours, France

* laurence.berville@imbe.fr or laurence-berville@hotmail.fr









7 workers from each nesting site (n =350) were immersed separately in 5µl of hexane (15 min) in order to extract and estimate cuticular hydrocarbons. 3µl of this extract was used for capillary gas chromatography (GC), carried out using a Varian 3 900 gas chromatograph equipped with a flame ionization detector (FID) and a Chrompack CPSil5WCOT apolar capillary column and interfaced with Star 5.5 (Varian) software. Oven temperature was held at 100°C, then increased to 220°C at 10°C/min, then to 320°C at 3°C/min and finally held at constant temperature for 10 min. The injector and flame-ionization detector were at 280 and 250 °C respectively. Quantitative data were obtained by integrating peaks.

Total body washes of 20 pooled workers from each species were immersed in 100µl of hexane (20 min) and used for extraction. 4µl of each extract were run into an Agilent 6890N GC equipped with Chrompack CPSil5WCOT apolar capillary column. The GC was coupled with a 5375 Agilent Technologies Mass Spectrometer (**GC-MS**). The carrier gas was helium at 1 ml/min with the injector in splitless mode. Oven temperature was isothermal at 70°C for 1 min, followed by 30°C/min to 180°C, then increased at a rate of 5°C/min to 320°C, and finally held for 15 min.


Results:






New PCA based on the proportion of the 77 CHCs of worker from the 3th group. K-means cluster analysis distinguishes 3 groups.

Cuticular profiles of the 5 Tapinoma species :





Quantitative presence in percent of the 5 distinct hydrocarbons groups according to the species.

T. s: 70 CHCs => C25 to TriMe C37

T. e: 39 CHCs => C25 to DiMe C31

T. is : 65 CHCs => C25 to TriMeC35

T. m : 62 CHCs => C25 to DiMe C35

 $T. N: 35 \text{ CHCs} \Rightarrow C27 \text{ to DiMe C32}$

Conclusions:

- The chemistry of the CHC bouquet is shown to be a good tool for taxonomists, being species-stable over thousands of kilometers.
- Extremely high diversity in the chemical recognition cues:
- Only 3CHCs in common to all species (C27, C29& 13 Me C31)
 - 108 CHCs (64.28%) present only in one species
- Only T. nigerrimum and T. israele present n-alkenes
- Across the 5 species, wide range of CHCs (168) were found, occuring between C25 to C 37



Max Planck Institute for Chemical Ecology

