

Structural and functional diversity of type IV secretion systems

Tiago R D Costa, Jonasz B Patkowski, Kévin Macé, Peter J. Christie, Gabriel Waksman

▶ To cite this version:

Tiago R D Costa, Jonasz B Patkowski, Kévin Macé, Peter J. Christie, Gabriel Waksman. Structural and functional diversity of type IV secretion systems. Nature Reviews Microbiology, 2024, 22 (3), pp.170-185. 10.1038/s41579-023-00974-3 . hal-04284264

HAL Id: hal-04284264 https://hal.science/hal-04284264

Submitted on 18 Dec 2023

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution - NonCommercial 4.0 International License

1	Structural and functional diversity of type IV secretion systems
2	
3	
4	Tiago R. D. Costa ^{1, *} , Jonasz B. Patkowski ¹ , Kévin Macé ^{2,3} , Peter J. Christie ^{4, *} and Gabriel Waksman ^{2, *}
5	
6	¹ Centre for Bacterial Resistance Biology, Department of Life Sciences, Imperial College, London, UK
7	² Institute of Structural and Molecular Biology, Birkbeck and UCL, London, UK
8	³ Institut de Génétique et Développement de Rennes (IGDR), Université de Rennes, CNRS, Rennes, France
9	⁴ Department of Microbiology and Molecular Genetics, McGovern Medical School at UTHealth, Houston, Texas,
10	USA
11	
12	
13	E-mails:
14	t.costa@imperial.ac.uk;
15	peter.j.christie@uth.tmc.edu;
16	g.waksman@ucl.ac.uk; g.waksman@bbk.ac.uk
17	
18	

19 Abstract

20 Considerable progress has been made in recent years in the structural and molecular biology of type IV secretion 21 systems (T4SSs) in Gram-negative bacteria. The latest advances have substantially improved our understanding 22 of the mechanisms underlying the recruitment and delivery of DNA and protein substrates to the extracellular 23 environment or target cells. In this Review, we aim to summarize these exciting structural and molecular biology 24 findings, and to discuss their functional implications for substrate recognition, recruitment and translocation as 25 well as the biogenesis of extracellular pili. We also describe adaptations necessary for deploying a breadth of 26 processes, such as bacterial survival, host-pathogen interactions, and biotic and abiotic adhesion. We highlight 27 the functional and structural diversity of this extremely versatile secretion superfamily to function under 28 different environmental conditions and in different bacterial species. Additionally, we emphasise the 29 importance of further understanding the mechanism of type IV secretion, which will support us in combating 30 antimicrobial resistance and treating T4SS-related infections.

31 [H1] Introduction

32 Type IV secretion systems (T4SSs) are a family of highly complex and versatile nanomachines that span the entire 33 cell envelopes of Gram-positive and Gram-negative bacteria as well as Archaea ¹⁻³ (Figure 1). They function in 34 two main capacities, as DNA transfer (conjugation) systems or as protein effector translocators^{4,5}, generally 35 mediating transfer of macromolecules by mechanisms requiring direct donor-target cell contact. A few systems 36 have evolved the capacity to export DNA or protein substrates to the extracellular milieu, or to take up DNA 37 from it. Many T4SSs also elaborate surface organelles such as conjugative pili or surface adhesins to promote 38 attachment and biofilm formation, but there are also examples of T4SSs that seem to have lost the capacity to 39 translocate substrates and instead function only in adhesion. Considering this enormous range of biological 40 activities, the T4SSs are exceptionally important from a medical perspective. Accordingly, they are increasingly 41 viewed as viable targets for therapeutic intervention to thwart the spread of conjugation-driven antibiotic 42 resistance and infection by pathogens ^{6,7}.

43

44 T4SSs in Gram-negative species are composed minimally of 12 core subunits that are generically termed VirB1 -45 VirB11 and VirD4 (Ref. 8). Systems assembled only with the core VirB–VirD4 components are considered 46 'minimized', and many of these systems function as conjugation machines by delivering DNA substrates to target 47 bacteria^{9,10}. Over the course of evolution, T4SSs have acquired several additional protein components that are 48 integrated into the core structure composed of VirB and VirD4 proteins. As a result, assembly of an expanded 49 T4SS may require up to 25 different proteins ^{10,11}. Some of these expanded systems can mediate conjugative 50 DNA transfer, but many have acquired new functionalities relating to translocation of effector proteins or toxins, 51 with or without retention of the ancestral DNA transfer function^{12,13}.

52

53 In recent years, our understanding of the architectures and mechanisms of actions of T4SSs has increased 54 substantially, most notably by implementation of state-of-the art microscopy techniques including cryo-electron 55 microscopy (cryo–EM), cryo-electron tomography (cryo–ET) and fluorescence microscopy. These studies have 56 enabled visualization of major machine subassemblies and conjugative pili at or near atomic resolution or fully 57 intact T4SSs in the native context of the bacterial cell envelope at lower resolutions (Supplementary Table 1). 58 Most exciting, a very recent study defined for the first time the architecture of a nearly completely intact 59 minimised T4SS encoded by the conjugative plasmid R388 (Ref. ¹⁴). There also has been considerable progress 60 in defining the architecture and mechanisms of action of the VirD4 components of T4SSs; these ATPases have 61 crucial roles in recruiting and coupling substrates to the translocation channel and hence are termed type IV 62 coupling proteins (T4CPs). Finally, studies applying a combination of *in situ* cryo-ET and fluorescence microscopy 63 have advanced our knowledge of T4SS assembly dynamics and spatial organization within intact cells.

64

This Review will primarily focus on the T4SSs found in Gram-negative bacteria, as they have been the most extensively characterized in terms of both structure and function. We will summarise the current knowledge of paradigmatic T4SSs functioning in these bacteria with a focus on their architectures and adaptations for specialised functions. We also update the reader on recent studies exploring the biogenesis pathways and spatial

3

69 localization of T4SSs, and we conclude with a brief review of progress toward developing small molecule 70

inhibitors of T4SSs and manipulating these versatile nanomachines for novel therapeutic ends.

71

72 [H1] Architectures of minimised systems

73 Early biochemical studies supplied evidence that the VirB subunits VirB7, VirB9 and VirB10 assemble as a 74 stabilizing structural scaffold for the T4SS; this scaffold ultimately was designated as the outer membrane core 75 complex (OMCC)^{8,15}. OMCCs are intrinsically stable and amenable to isolation for structural characterization. 76 Accordingly, over a decade ago, high-resolution structures were presented for the OMCC associated with a 77 minimized T4SS elaborated by the conjugative plasmid pKM101 (T4SS_{pKM101})^{16,17}, and soon afterwards for several 78 other OMCCs from minimized systems¹⁸⁻²¹. The most recent structure presented for the nearly intact T4SS 79 encoded by plasmid R388 (T4SS_{R388}) now has provided important refinements of these earlier structures ¹⁴ 80 (Figure 2a). The OMCC from the T4SSR388 presents as a barrel-shaped structure of 130 Å in height and 185 Å in 81 width. It is composed of an outer and inner layer that are designated as the O-layer and I-layer, respectively,. 82 Remarkably, in contrast to the lower-resolution structures obtained previously, the new structure shows that 83 the O- and I-layers have different symmetries (Figure 2a). The O-layer is made up of 14 copies of homologs of 84 the VirB7 lipoprotein and C-terminal domains (CTDs) of the VirB9- and VirB10-like subunits. VirB7 anchors the 85 OMCC to the outer membrane via its N-terminal lipid modification. The I-layer assumes a 16-fold symmetry, 86 formed by copies of N-terminal domains (NTDs) of the VirB9 and VirB10 subunits. The C14:C16 symmetry 87 mismatch between the O- and I-layers is accommodated by a unique configuration of two VirB9 and VirB10 88 subunits; while their NTDs insert in the I-layer, their CTDs do not form part of the O-layer ¹⁴ (Figure 2a). 89 Importantly, the VirB10_{CTD} of the O-layer contains a hydrophobic 2-helix bundle (termed the antennae 90 projection) extending from the VirB10_{CTD} β-barrel domains, and 14 of these bundles assemble to form a hollow 91 pore through the outer membrane. This pore has a diameter of 32 Å at the extracellular entrance, which is 92 proposed to expand further through a hinge-like conformational change between the antennae projection 93 helices¹⁴ to accommodate a growing pilus and substrate translocation. This mode of action is in line with earlier 94 findings that VirB10 associated with the Agrobacterium tumefaciens VirB-VirD4 T4SS (T4SSAgro) undergoes a 95 conformational switch in response to sensing and transduction of intracellular signals to gate the outer 96 membrane channel ²²⁻²⁵.

97

98 The phytopathogen Xanthomonas citri elaborates a T4SS whose constituent subunits are close homologs of 99 components of the pKM101- and R388-encoded conjugation machines, but this system (T4SS_{x. citri}) functions as 100 an interbacterial killing machine by delivering toxins to neighbouring bacteria²¹. Interestingly, the OMCC of the 101 T4SS_{X. citri} does not display a mismatch between the O- and I- layers as both assume a 14-fold symmetry ²¹. It 102 encodes a much larger VirB7 homologue with a pronounced C-terminal N0 domain ²⁶, a feature also found in 103 the Dot/Icm system (T4SS_{Dot/Icm}) encoded by Legionella pneumophila (see below)^{27,28}. The NO domain widens 104 horizontally, giving rise to a flying saucer-shaped structure rather than the barrel-like structures of OMCCs from 105 the T4SS_{pKM101} and T4SS_{R388} machines. Whether or how these structural variations contribute to the specialized 106 function of the T4SS_{X, citri} machine in interbacterial transmission of toxins is not yet known.

107

108 The nearly intact T4SS_{R388} also details the organization of the inner membrane complex (IMC) and periplasmic 109 stalk (Figure 2a). The stalk connecting the OMCC and IMC has an overall diameter of ~92 Å and length of 216 Å. 110 The stalk consists of a pentamer of VirB5 subunits located proximally to the OMCC connected to a pentamer of 111 VirB6 that inserts into the inner membrane by an N-proximal α -helix ¹⁴ (Figure 2a). Six knobs (or 'arches') 112 composed of the C-terminal domains of VirB8 subunits surround the central stalk near the inner membrane. The 113 architecture of the stalk is intriguing in view of prior evidence that the VirB5 subunit localizes at the tip of the T-114 pilus elaborated by the A. tumefaciens VirB/VirD4 system ²⁹. Moreover, the finding that the bulk of VirB6 resides 115 in the periplasm completely revises our views of how VirB6 contributes to substrate translocation and pilus 116 assembly. VirB6 subunits are highly hydrophobic, and previously were envisioned to adopt a polytopic topology 117 in the inner membrane and form part of the inner membrane channel ^{1,30}. In the T4SS_{R388} structure, the bulk of 118 VirB6 is shown to assemble as a pentameric platform entirely in the periplasm while only one N-terminal 119 hydrophobic domain made of two trans-membrane helices integrates in the inner membrane (Figure 2a). Thus, 120 most of the hydrophobic domains of VirB6 subunits are shielded from the aqueous environment of the periplasm 121 through extensive intersubunit contacts. It is also noteworthy that the pentameric symmetry of the VirB5–VirB6 122 stalk matches that of the conjugative pilus.

123

124 Using structural-complementarity simulations, the VirB2 pilin component of the conjugative pilus was 125 successfully docked on the VirB6 platform, which led to a new model for how conjugative pili assemble (Figure 126 2b). Upon synthesis, pilin subunits integrate into the inner membrane, forming a pool for subsequent use in 127 building the pilus³¹ (Figure 2b, step 1). Co-evolution studies confirmed by site-directed mutagenesis identified 128 one of the two trans-membrane helices as binding and recruitment site for the VirB2 pilin. Upon receipt of an 129 unknown signal, five pilin subunits are extracted from the inner membrane into the periplasmic assembly site 130 located in the interacting surfaces of the VirB6 and VirB5 proteins (Figure 2b, step 2). Through reiterative 131 recruitment and extraction of pilins to the VirB6 platform, the pilus extends and displaces the VirB5 pentamer 132 upwards and through the outer membrane to the extracellular milieu ¹⁴ (Figure 2b, step 3). The VirB5 pentamer 133 remains at the tip of the pilus, where it binds specific receptors on the target cell surface or embeds directly into 134 the target cell membrane ³² (Figure 2b, step 4). The role of VirB5 as an adhesin was initially speculated in 2008 135 (Ref. ³¹), and subsequent evidence has been found in the Cag system. In this system, the functional orthologue 136 of VirB5, termed CagL³³⁻³⁵, has been shown to bind integrin on the human cell surface ³². Furthermore, VirB5 has 137 exhibited conformational changes similar to those observed in pore-forming proteins ¹³, suggesting potential 138 structural parallels with the translocon pores of type III secretion systems.

139

The IMC of the T4SS_{R388} is dominated by two concentric rings that extend into the cytoplasm (Figure 2a). These rings comprise a hexamer of dimers of the VirB4 ATPase ¹⁴. The dimers are arranged so that one protomer forms the inner ring and the other forms the outer ring. Protomers of each dimer are connected by their NTDs, which also embed into the inner membrane. Additionally, the NTDs of protomers comprising the inner and outer hexameric rings form specific contacts with the VirB3 and VirB8 subunits, respectively. This interaction network

anchors the VirB4 hexamer of dimers in the inner membrane and potentially enables ATP-dependent structural
 changes ¹⁴ (Figure 2a).

147

148 Structures of the T4SS_{pKM101} and of several mutant machines assembled within the cell envelope also were solved 149 by in situ cryo-ET ³⁶ (Figure 2a). Despite their lower resolutions, these structures revealed several distinctive 150 features. For example, during assembly of the T4SS_{pKM101}, the outer membrane is extensively remodelled as 151 evidenced by invagination of the outer leaflet and the absence of the inner leaflet of the OM at the machine -152 OM junction. The OMCC possesses a central chamber sufficiently large to accommodate a growing pilus as it 153 extends from the stalk assembly platform. Most importantly, analyses of the IMCs from WT and mutant strains 154 established that the VirB4 subunit assembles in vivo as a central hexamer of dimers at the channel entrance, 155 and that neither of the VirB11 or VirD4 ATPases contribute to IMC densities ³⁶. These latter ATPases also were 156 not visualized in the in vitro T4SS_{R388} structure, but the VirB11–VirB4 interaction was characterised using 157 computational methods and validated biochemically and by site-directed mutagenesis. VirB11 and VirD4 might 158 dock transiently with the VirB channel in response to unknown signals ¹³.

159

160

161 [H1] Architectures of expanded systems

Expanded T4SSs are composed of up to 25 distinct proteins, encompassing the 12 core components of VirB and VirD4 along with additional protein subunits ^{10,11}. The most widely studied expanded systems are the F plasmidencoded conjugative machine, and the *L. pneumophila* Dot/Icm and *Helicobacter pylori* Cag effector translocator systems. Although detailed structures for the IMCs from these systems remain to be determined, near atomic resolution structures are now available for OMCCs from all three systems, F plasmid-encoded conjugative pili, and the T4SS_{Dot/Icm} substrate recognition platform.

168

169 The OMCC associated with the F system is composed of the three core VirB components, VirB7-like TraV, VirB9-170 like TraK and VirB10-like TraB ^{15,37} (Figures 3a and 3b). It adopts a 'flying saucer' shape and is much larger than 171 the OMCCs of the minimized conjugation machines, with a diameter of ~268 Å and width of 115 Å. The complex 172 consists of concentric rings termed the outer and inner rings. The inner ring exhibits a 17-fold symmetry and is 173 composed of seventeen CTDs of TraB and NTDs of TraV. The outer ring displays a 13-fold symmetry and is 174 composed of 26 CTDs of TraV and 26copies of TraK (Figure 3a). The TraK subunits assemble as 13 dimers, with 175 the CTDs of each dimer pair forming 13 elongated knobs that extend radially from the centre of the complex. 176 The TraB and TraV proteins were observed to form flexible linkers connecting the symmetrically mismatched 177 inner ring and outer ring. Structural flexibility imparted by the symmetry mismatch could account for the 178 dynamic properties of the F pilus during extension and retraction. This ability could also be potentiated by a 179 number of additional proteins present in the expanded F system compared to minimised conjugative systems, 180 which are not capable of pilus retraction.

181

182 The L. pneumophila Dot/Icm and H. pylori Cag systems function as effector translocators that aid in infection 183 processes (Figure 1). The Dot/Icm system has the remarkable capacity to translocate at least 330 effector 184 proteins, nearly 10 % of the proteome, into eukaryotic cells ^{28,38,39}. Effector translocation induces a myriad of 185 physiological changes marked by conversion of human phagosomes into replication-permissive compartments 186 called Legionella-containing vacuoles (LCV)⁴⁰. The OMCC_{Dot/Icm} has a highly complex organisation and is 187 considerably larger (~420 Å in diameter) than that of the F system (Figure 3a). VirB7-like DotD, VirB9-like DotH, 188 and VirB10-like DotG proteins form the structural scaffold, which is built upon by incorporation of other Dot/Icm-189 specific proteins including DotF, DotC, Dis1, Dis2 and Dis3 (Refs. ^{41,42}). The OMCC_{Dot/Icm} has a three-layer topology 190 with a disk-shaped outer membrane cap (OMC) of 13-fold symmetry, a proximal dome of 16-fold symmetry and 191 a smaller periplasmic ring of 18-fold symmetry (Figure 3a). Distinct O- and I- layers reminiscent of the OMCC_F 192 are likely to be represented by the OMC and PR domains, respectively. In line with this proposal, VirB9-like DotH 193 was observed to comprise both the OMC and periplasmic ring reminiscent of the architectures of the VirB9 194 counterparts in the F, pKM101 and R388 systems.

195

196 The Cag system of *H. pylori* secretes the oncoprotein CagA as well as several other nonproteinaceous substrates 197 into human epithelial cells ⁴³⁻⁴⁵ (Figure 1). These T4SS_{Cag}-mediated functions induce pathological changes in 198 epithelial cells that facilitate H. pylori infection of the gastrointestinal tract ⁴⁶. The T4SS_{Dot/Icm} and T4SS_{Cag} are the 199 largest of the T4SSs characterized to date and both possesses some of the structural properties identified in the 200 other systems. The mushroom-shaped OMCC_{cag} is 400 Å in width and 250 Å in height, and is composed VirB7-201 like CagT, VirB9-like CagX and VirB10-like CagY plus two system-specific proteins, Cag3 and CagM ⁴⁷ (Figure 3a). 202 The large size of the OMCC_{Cag} is conveyed by the larger sizes of CagT, CagX and CagY relative to their VirB 203 counterparts, and by incorporation of multiple copies of Cag3 and CagM in the structure. The OMCC_{Cag} consists 204 of an OMC with clear O- and I-layers. The OMC is made up of CagT, CagX, CagY, Cag3 and CagM, and the 205 periplasmic ring composed of CagY and CagX. Portions of the periplasmic ring show structural resemblance to 206 the I-layer of OMCC_{X. citri} (Ref. ⁴⁸). Both the inner and outer regions of the OMC have a 14-fold symmetry, whereas 207 the periplasmic ring assumes a 17-fold symmetry (Figure 3a). As shown for other VirB9 subunits, CagX forms 208 contacts with both the OMC and periplasmic ring, thus bridging the symmetry mismatch between those two 209 complexes ⁴⁸. Reminiscent of the F system, the observed asymmetries among OMCC substructures of the L. 210 pneumophila Dot/Icm and H. pylori Cag systems might have evolved to confer specialized functions, in these 211 cases relating specifically to the infection process. OMCC structural flexibility might, for example, be needed to 212 establish dynamic yet productive contacts with eukaryotic host cells or to coordinate the timing and delivery of 213 substrates into the eukaryotic cell targets.

214

These expanded T4SSs also have been visualized by *in situ* cryo-ET in their native cellular contexts. Remarkably, F systems elaborate a presumptively quiescent translocation channel and three morphologically distinct platforms upon which the F pilus is docked⁴⁹ (Figure 3b). The channel, designated as the F1-channel complex, consists of the OMCC joined to the IMC by a periplasmic cylinder that is thicker and more pronounced than the thin periplasmic stalk of the T4SS_{R388} (Figure 3b). The central hexamer of dimer configuration of the VirB4-like

220 TraC ATPase is readily visualized at the cytoplasmic entrance, reminiscent of the in situ T4SS_{DKM101} and recent in 221 vitro high resolution T4SS_{R388} structures (Figure 2a). An F-pilus associated structure, termed the F2-channel-222 pilus complex, resembles the F1 complex but has the F pilus attached at the distal end of the OMCC (Figure 3b). 223 The F1 and F2 complexes are postulated to correspond to the quiescent channel and active pilus-assembly 224 factory involved in mate-seeking and mating. Two other F pilus-associated structures, designated as the F3-talk-225 pilus and F4–outer membrane-pilus complexes, consist of the F pilus attached respectively to a thin stalk density 226 that spans the periplasm or a small outer membrane density without any underlying structure. The F3 and F4 227 complexes lack discernible channels for substrate transfer or pilus assembly, and accordingly are proposed to 228 function exclusively as holding platforms for nonretractile F pili. These inert structures might contribute to 229 nonspecific cell aggregation and biofilm formation or as decoys for bacteriophages that rely on F pilus retraction 230 to gain access to the cell envelope 49.

231

232 Visualization of the in situ T4SS_{Dot/Icm} structure revealed an OMCC with a 'Wi-Fi symbol'-like architecture of 400 233 Å in diameter whose assembly is dependent on the DotC, DotD, DotF, DotG and DotH subunits ^{50,51} (Figure 3c). 234 The OMCC is connected to a cylinder with an outer diameter of ~20nm and a central lumen or channel of ~6 nm, 235 which narrows to a diameter of ~10nm and a channel of ~3nm near the IM⁵¹. As with other T4SSs (Refs. ^{14,21,36}), 236 at the cytoplasmic face of the IMC is the hexamer of dimer configuration of VirB4 ATPase (DotO). Remarkably, 237 VirB11-like DotB, a second ATPase of this system, assembles as a hexamer that dynamically associates with the 238 DotO inner hexamer by a mechanism dependent on cycles of ATP binding and hydrolysis (see below)⁵¹. Recently, 239 another Dot/Icm system elaborated by Coxiella burnetti has been determined by a combination of cryo-ET and 240 cryo-focused ion beam (cryo-FIB) milling to closely resemble that of L. pneumophila 52. Interestingly, this study 241 documented a correlation between assembly of the machinery and developmental transitions of C. burnettii 242 cells during infection, which complements previous observations that reported a dependence of effector 243 translocation by this T4SS on the progression of *C. burnetti* infection⁵³.

244

245 In the *in situ* cryo-ET map of the *H. pylori* Cag machine, the OMCC resembles the equivalent substructure solved 246 at a higher resolution by cryo-EM ^{48,54-56}. Analyses of mutant machines confirmed that in situ assembly of the 247 OMCC requires the Cag subunits CagX, CagY and CagM, and that CagT and Cag3 contribute to peripheral 248 densities. Other noteworthy features of the in situ T4SS_{Cag} machine include a hollow cylinder that extends across 249 the periplasm, connecting the OMCC to the IMC. The IMC is architecturally more complex than the equivalent 250 substructures of the F and Dot/Icm systems in having three concentric rings instead of two. A unique feature 251 among T4SSs is that the extracellular domain of the CagY subunit of the T4SS_{Cag} contains multiple binding sites 252 for Toll-like receptor 5 (TLR5) and functions in regulating immune responses of the host⁵⁷. The inner and middle 253 rings correspond to the central and outer hexamers of VirB4-like CagE, but assembly of the outer ring is 254 dependent on production of VirD4-like Cag8⁵⁵. VirB11-like Caga associates at the base of the CagE central 255 hexamer, but in contrast to the Dot/Icm system, Caga seems to associate stably with CagE and not dynamically 256 as a function of ATP binding and hydrolysis as shown for DotB⁵⁵ (Figure 3d). A comprehensive table that

8

257

summarises all available structures to date comprising both, the minimised and expanded systems is provided as Supplementary Table 1.

258 259

260

261 [H1] Other T4SS machine adaptations

262 Besides appropriating novel subunits for functional diversification, many T4SSs of both the minimised and 263 expanded types have diversified through modifications of core VirB components. VirB2-like pilins or pili or 264 associated surface adhesins have been adapted to enhance T4SS targeting to specific cell types (see below), 265 whereas certain IMC or OMCC components have undergone modifications through acquisition of novel domains 266 for broadened T4SS functionality (Figure 4). These modifications generally function to direct machine assembly 267 at specific sites within the cell or coordinate T4SS localization or function with the cell cycle, or to promote 268 attachment to specific target cells, as described in more detail in supplementary text box 1. Symmetry 269 mismatches are likely to provide regions of mobility or flexibility between different protein complex layers or 270 subassemblies. However, despite the presence of symmetry mismatch in several types of secretion systems, 271 precisely how symmetry mismatch contributes to machine functions at mechanistic or structural levels have not 272 been established for any characterized secretion system ¹⁰.

273

274 The VirB6 subunits are the most extensively modified of the IMC components. In many T4SSs, VirB6 subunits are 275 composed predominantly of 5-7 hydrophobic domains that are likely to assemble as part of the central stalk 276 structure, as demonstrated for the R388-encoded VirB6 subunit (Figure 1a). However, many VirB6 subunits 277 termed 'extended-VirB6s' are considerably larger by virtue of the presence of one or more large central of C-278 terminal hydrophilic domains¹ (Figure 4). Remarkably, where characterized these hydrophilic domains 279 contribute in distinct ways to establishment or inhibition of productive donor-target cell interactions. In F 280 systems, a large ~600 residue C-terminal hydrophilic domain (CTD) of VirB6-like TraG is involved in entry 281 exclusion, a process that blocks redundant DNA transfer between the donor cells. When donor cells form mating 282 junctions with other donor cells, the CTD of TraG produced by one donor cell establishes contact with TraS_F, an 283 inner membrane protein produced by the paired donor cell ^{58,59}. This contact is achieved either by extension of 284 TraG_F across the mating junction or by proteolytic cleavage and translocation of the CTD of TraG_F through the 285 T4SS_F into the paired donor cell. In *Rickettsia spp.*, multiple copies of extended-VirB6s are present with sizes 286 ranging from 600 to over 1400 residues ^{60,61}. Large hydrophilic domains are surface displayed where they are 287 implicated in promoting endosymbiotic or pathogenic relationships with eukaryotic target cells⁶². In the L. 288 pneumophila Dot/Icm system, VirB6-like DotA possesses multiple hydrophobic domains flanking a central 289 hydrophilic domain. Although DotA associates with the inner membrane reminiscent of other VirB6 family 290 members, it can also be exported in a Dot/Icm T4SS-dependent manner to the extracellular milieu where it forms 291 ring-like oligomers whose functions are presently unknown ⁶³.

292

All three of the core VirB components of OMCCs can be modified through acquisition of novel motifs (Figure 4).
 Although many VirB7-like lipoproteins are small (~50 residues) and resemble the archetypal *A. tumefaciens*

VirB7, many others are considerably larger (~160-300 residues) as shown for *H. pylori* CagT, *X. citri* VirB7 and *L. pneumophila* DotD ^{11,21,41,48}. In the *H. pylori* system, both CagX and CagY are considerably larger than their VirB9 and VirB10 counterparts. Indeed, in the case of CagY, only the extreme C-terminal region adopts the characteristic VirB10 β-barrel-domain folds that assemble as the central rings of OMCCs. A large middle repeat region (MRR) is composed of multiple repeats, which undergo extensive rearrangement during gene expression. These rearrangements yield many CagY variants that have been shown to regulate T4SS_{Cag} function — positively or negatively — to maximize persistent infection ⁶⁴⁻⁶⁶.

302

303

304 [H1] VirD4 substrate recruitment and translocation

For most T4SSs, VirD4-like ATPases also known as type IV coupling proteins (T4CPs) ^{22,23,67} are responsible for 305 306 substrate recruitment. T4CPs characteristically possess an N-terminal transmembrane domain (TMD) implicated 307 in establishment of contacts with IMC components of cognate T4SSs. A conserved nucleotide-binding domain 308 (NBD) is thought to provide the energy for early-stage substrate processing; for example, unfolding and opening 309 of the channel for substrate transfer. Two sequence-variable domains, the all-alpha domain (AAD) and C-310 terminal domain (CTD) if present, contribute to substrate recruitment. An X-ray structure of R388-encoded 311 TrwB, currently the structural archetype for T4CPs, showed that the NBD assembles as a homohexamer ⁶⁸ (Figure 312 5a). This NBD architecture bears similarities to the FtsK and SpolllE families of DNA motor proteins involved in 313 DNA translocation during cytokinesis ⁶⁷. The AAD sits at the base of the hexamer, optimally positioned for 314 docking of secretion substrates. Indeed, evidence has now been presented for binding of AADs to secretion 315 substrates in the A. tumefaciens and Xanthomonas VirB-VirD4 systems^{69,70}. Protein substrates are recognised 316 by the T4SS through conserved motifs called translocation signals. These highly specific signals encoded by the 317 protein substrates are tailored for the timely engagement and controlled secretion by the T4S apparatus. Details 318 on the different translocation signals associated with T4SS substrates can be found in the supplementary text 319 box 2.

320

321 VirD4 homologues also interact with VirB11 ATPases ^{22,71,72}. VirB11 family members are structurally similar 322 (Figure 5b) insofar as their NTDs and CTDs resemble each other and are connected by a flexible linker. VirB11 323 subunits assemble as homohexamers⁷³, and in this oligomeric state, the interdomain linkers facilitate fluent 324 domain swaps without affecting hexamer assembly⁷⁴. The VirB11, VirD4 and VirB4 ATPases act in concert to 325 orchestrate pilus biogenesis and substrate transfer, although mechanistic details underlying this coordination of 326 ATPase functions remain unknown ^{22,71,72,75,76}. The VirB11 ATPases are viable drug targets, as illustrated by an 327 early report that small molecule inhibitors of VirB11-like Cagα blocks virulence of *H.pylori*⁷⁷.

328

329 In the *L. pneumophila* Dot/Icm system, the CTD of VirD4-like DotL has a highly complex role in recruitment of its 330 many protein effectors and their subsequent secretion through the T4SS_{Dot/Icm} apparatus (Figure 5c). Although 331 recruitment of various effectors to DotL can proceed independently of known associated chaperones, 332 recruitment of others is strictly dependent on the IcmS and IcmW chaperones, with or without an additional

10

333 requirement for the LvgA chaperone⁷⁸⁻⁸¹. DotL assembles as a hexamer, and associated with each of the six CTDs 334 are various Dot adaptors including DotM, DotN, DotY, DotZ, IcmS and IcmW. The DotM adaptor itself also acts 335 as a recruitment platform for the some of the IcmSW-independent effectors⁸². The six CTD-adaptor complexes 336 together form a bell-shaped structure, termed the type IV coupling complex (T4CC), that extends into the 337 cytoplasm and functions as the substrate recruitment platform (Figure 5c). The T4CC possesses at least two 338 known effector binding-sites, one on DotM and a second strictly dependent on the IcmS and IcmW chaperones, 339 with or without an additional requirement for the LvgA chaperone⁷⁸⁻⁸¹. Through binding of distinct arrays of 340 effectors based on their associations with different chaperones, the T4CC is envisioned to regulate substrate 341 transfer during *L. pneumophila* infection⁸³⁻⁸⁸.

342

343 In the *H. pylori* Cag T4SS, VirD4-like Cagß is structurally similar to TrwB (Figure 5a); however its AAD contributes 344 in a unique way to modulation of Cag β function. In this system, two cytosolic proteins, CagF and CagZ, function 345 together with Cag β to recruit the CagA substrate to the T4SS_{Cag}. Prior work showed that CagF functions as a 346 chaperone by binding a 100-residue region of CagA, whereas CagZ stabilizes Cag β ⁸⁹. Recently, crystal structures 347 presented for the Cag β AAD–CagZ interaction led to a model whereby the CagZ–AAD contact maintains Cag β in 348 a monomeric state, thereby suppressing ATPase activity and rendering Cag β inactive⁹⁰. Upon receipt of an 349 unknown signal, CagF binds and recruits CagA to Cag β , resulting in release of CagZ and assembly of the Cag β 350 hexamer. This catalytically active form of Cag³ then unfolds and translocates CagA through the T4SS_{Cae}.

351

352 Among the conjugative T4SSs, T4CPs recruit a specialised DNA-processing complex called the relaxosome ^{5,91} 353 (Figure 6a). In the well-characterized F system, the relaxosome is composed of four proteins assembled at the 354 origin-of-transfer (oriT) sequence of the F plasmid. The largest relaxosome subunit, Tral relaxase (~200 kDa), 355 possesses a trans-esterase domain that nicks and covalently attaches to the 5' end of the nicked DNA strand 356 destined for translocation (T-strand). Tral also possesses a vestigial helicase domain that serves as an single-357 stranded DNA (ssDNA)-binding domain, an active 5' to 3' helicase domain, and a C-terminal domain that 358 functions as a recruitment platform for the remaining relaxosome components^{92,93}. Two other relaxosome 359 components TraY and IHF cause conformational changes in the DNA topology that expose the nick (nic) site for 360 cleavage by Tral. ^{94,95}. Finally, TraM is a homotetramer responsible for docking the DNA–relaxosome complex to 361 the coupling protein TraD ^{96,97} (Figure 6a). On the DNA side, a pair of TraM tetramers recognise a particular DNA 362 sequence known as *sbmABC* motifs within the *oriT* of the F-plasmid ⁹⁸, resulting in up to six tetramers bound to 363 one F plasmid. Each TraM tetramer is subsequently recognised by C-terminal tails of the hexameric TraD⁹⁶ (Figure 364 6a). This interaction ensures F plasmid docking and transfer while actively blocking translocation of a co-existing 365 conjugative plasmid(s) through the F machine ^{84,97,98}.

366

Remarkably, at this time there is little detailed structural information concerning the physical relationship
 between the T4CP and the T4SS channel. Consequently, the route or routes by which DNA and protein substrates
 are conveyed through T4SSs is not known, and two distinct translocation pathways have been envisioned (Figure
 As illustrated for the T4SS_{Dot/Icm} machine, upon recruitment of substrates to the DotL T4CC, the T4CC functions

371 in one of two ways, which ultimately dictates whether substrates are delivered in one or two steps across the 372 entire cell envelope (Figure 6B). In the one-step translocation pathway⁹⁹, DotL is situated beneath or near the 373 T4SS and, upon substrate engagement, the T4CC passes captured substrates to the base of the T4SS channel 374 marked by the DotB and DotO hexameric ATPases. The DotL-DotB-DotO ATPase ternary complex then 375 orchestrates substrate unfolding, dissociation of chaperones and/or adaptors, and delivery of the translocation 376 intermediate into the channel for conveyance in one step through a channel that extends from the cytoplasmic 377 face of the inner membrane to the cell exterior (Figure 6b). In the alternative two-step translocation model ⁹⁹, 378 DotL is situated in physical proximity to the T4SS. The T4CC captures and then shunts substrates directly into the 379 lumen of the DotL hexamer (whose NTD spans the inner membrane) for delivery across the inner membrane. 380 Once in the periplasm, in a second translocation reaction, substrates are recruited to and enter the T4SS channel 381 for delivery to the cell surface ^{83,51,100} (Figure 6b). For reasons outlined in the next paragraph, we hypothesize 382 that delivery of the F plasmid transfer intermediate follows the one-step pathway, whereby TraD docks at the 383 base of the channel in complex with VirB4-like TraC⁸⁴ (Figure 6a).

384

385 In the absence of detailed structural information about the T4CP–T4SS connection, we currently favor the one-386 step pathway for DNA as well as protein substrates, given the potentially deleterious consequences 387 (degradation, misfolding, misrouting and temporal disruption) encountered by ssDNA transfer intermediates or 388 as many as several hundred effectors delivered into the bacterial periplasm via the two-step pathway. Early 389 crosslinking studies in A. tumefaciens system also supplied experimental support for a one-step pathway. In 390 those studies, DNA substrates of the VirB–VirD4 system were shown to engage in sequential order with ATPase 391 subunits comprising the cytoplasm-inner membrane interface, then with components of the IMC, and finally with components of the OMCC and pilus ^{22,23,30,71,101}. Structural studies of the T4SS_{R388} by negative stain-EM ¹⁰² 392 393 and of the T4SS_{Cag} by *in situ* cryo-ET⁴⁹ also have supplied evidence that T4CPs can associate with the VirB4 394 hexameric platform at the channel entrance, thus being optimally positioned to deliver captured substrates 395 directly into the channel. While there are clear examples of T4SSs that utilize two-step translocation pathways 396 to deliver substrates to the extracellular milieu ¹⁰³ or to eukaryotic target cells ¹⁰⁴, these systems lack VirD4 T4CPs 397 and thus rely on alternative secretion systems for recruitment and translocation of substrates across the inner 398 membrane.

399

400

401 [H1] Conjugative pili and target cell attachment

T4SSs elaborate various surface structures that have important roles in promoting donor-target cell contacts.
Among these, conjugative pili are major contributors to the rapid and widespread dissemination of plasmids and
other mobile elements, and their cargoes of antibiotic resistance determinants, among Gram-negative bacteria.
High-resolution structures recently have been generated for several conjugative pili, including those produced
by F plasmids (F pilus)^{105,106}, the *A. tumefaciens* VirB-VirD4 system (T pilus)¹⁰⁷⁻¹⁰⁹ and IncN plasmid pKM101¹⁰⁸.
Throughout evolution, the relationships between hosts and pathogens have led to adaptations in T4SSassociated pili, resulting in more specialized functions. Using a 'one-size-fits-all' approach is not effective for

12

adhesion to distinct host-encoded receptors or different biotic and abiotic surfaces. As a result, pili and pilins
have diversified to fulfil these functions. The supplementary text box 3 summarizes the diversification of pili and
bacterial adhesins and their role in establishing close contacts with target cells.

412

Structurally, assembled pili differ in rotational raise of subunits (Figure 7a) and can display either a five-start helical symmetry, as shown for T- and N-pili and F-pili encoded by the classical F plasmid and F-like pED208, or a one-start symmetry as shown for another F-pilus encoded by pKpQIL or an Archaeal conjugative pilus¹⁰⁹. The outer diameters of conjugative pili in Gram-negative species are 76 Å - 87 Å and inner lumens are 23 Å - 26 Å in width, whereas the archaeal pilus has an inner diameter of 16 Å (Figure 7b).

418

419 As noted earlier, pilin subunits accumulate in the inner membrane as a pool for recruitment to build the 420 conjugative pilus upon receipt of an unknown signal. Intriguingly, structural studies of the conjugative pili have 421 established that, during extraction from the inner membrane, pilins co-extract phospholipid molecules. In F 422 systems, for example, TraA pilins co-extract phosphatidylglycerol molecules in a 1:1 stoichiometry (Figure 7C), 423 resulting in assembly of the helical fiber with phosphatidylglycerol molecules lining the pilus lumen. These 424 phosphatidylglycerol molecules impart an overall negative charge to the F pilus lumens^{105,106}. The A. tumefaciens 425 T pilus^{108,109} and pKM101-encoded N pilus ¹⁰⁸ share the general features first reported for F pili, but the lumens 426 are lined with phospholipids with different head groups (Figure 7c). For N pili this results in an overall negative 427 charge reminiscent of F pili¹⁰⁸, but for T pili the lumen has an overall positive charge. Another interesting 428 difference is the presence of a kink between the first alpha helices of pilins comprising the T- and N- pili, which 429 is absent in pilins associated with pili from expanded systems (Figure 7c; arrow). Why conjugative pili differ in 430 phospholipid compositions and what is the significance of the kink between helices $\alpha 1$ and $\alpha 2$ is unknown, but 431 it is reasonable to propose that these features impart biophysical properties of importance for specialized 432 functions in different environmental conditions.

433

434 In support of this model, recent biophysical studies have demonstrated that F pili are extremely flexible and 435 have spring-like properties with pronounced structural and thermochemical robustness. These properties are 436 postulated to accelerate conjugation rates and biofilm formation by F-plasmid-carrying cells by allowing for 437 effective function even in highly turbulent environments, such as those present in human gastrointestinal tracts 438 ¹¹⁰. This unique feature of F-pili might well be responsible for recent evidence that IncF plasmids are the most 439 dominant types of conjugative plasmids present in enterobacterial isolates from humans and animals ¹¹¹. In 440 contrast to F pili, other conjugative pili are generally shorter and more rigid and have not been shown to 441 retract¹¹². Because of their distinctive biophysical properties, these rigid pili are envisioned to readily break from 442 the cell surface, accumulate in the milieu, and mediate nonspecific aggregation of donor and potential recipient 443 cells, thus acting indirectly to facilitate propagation of mobile elements.

- 444
- 445

446 [H1] Conclusions and outlook

447 The remarkable recent progress in defining the architectures of T4SSs, both *in vitro* and *in situ*, confirms that the 448 functional diversity of these fascinating nanomachines is recapitulated at the structural level. Although there 449 clearly are emergent structural themes, most if not all T4SSs also have acquired systems-specific properties. This 450 is particularly evident with the expanded systems, which have appropriated components from unknown 451 ancestries that physically enlarge and add structural complexity to the T4SSs. The contribution of T4SS structural 452 diversity to their varied functions is still not fully understood. However, it is possible to speculate on potential 453 relationships. For instance, differences in pilin or adhesin subunits may result in variations in pilus biogenesis, 454 assembly and adherence properties. Another factor to consider is cargo recruitment, where structural 455 differences in the T4CP receptor or complex are likely to have evolved to specify restricted or expanded 456 substrate repertoires. Additionally, T4SS diversity at both the structural and functional levels are likely to have 457 evolved for recognition of specific host cell receptors, ultimately dictating the range of hosts that the T4SS can 458 target. Interestingly, DNA and protein substrates also have evolved a bewildering array of translocation signals 459 (supplementary text box 2), along with deployment of chaperones or adaptors, for docking with VirD4 substrate 460 recruitment platforms. In conjunction, VirD4 evolved to carry sequence-variable motifs or domains (for example, 461 AAD, CTD and T4CC) to specify and temporally regulate loading and delivery of substrates into the translocation 462 channel.

463

464 The T4SS field is well-poised to answer several long-standing questions. Most importantly, we still lack critical 465 information about the physical and functional relationship of the VirD4 T4CP with the cognate T4SS channel. 466 Also, what constitutes the translocation channel across the periplasm and the route of substrate transfer across 467 the entire cell envelope? At the cell exterior, how do mating junctions form, how are they physically configured, 468 and how do they dissociate after substrate transfer is completed? To what extent or for which systems do 469 conjugative pili routinely transmit DNA and protein substrates through their lumens? To address these 470 questions, a combination of structural analyses of mutant machines, fixed in their activated states or actively 471 engaged in substrate transfer, and other approaches such as correlative light and electron microscopy (CLEM) 472 offer considerable promise. Finally, and most central to this Review, how do the various T4SS structural 473 adaptations confer system-specific functions? A full answer to this question will be generated by widening the 474 current subjects of study to include the many other T4SSs functioning in diverse species of bacteria and archaea. 475

476 There is growing interest in translational initiatives aimed at blocking or repurposing T4SSs for therapeutic ends. 477 High-throughput screens are being used to identify small-molecule inhibitors, with goals of blocking conjugative 478 dissemination of antibiotic resistance or inhibiting effector translocation to suppress pathogenesis^{6,7,113-120}. 479 Structural advances of the different T4SSs summarized in this Review will continue to facilitate the rational 480 design of small molecules effective at blocking critical subunit–subunit interfaces. Conversely, as T4SSs are the 481 only bacterial secretion systems capable of delivering DNA or protein substrates to a wide range of target cell 482 types, these nanomachines are excellent delivery systems for therapeutic interventions. Indeed, the early 483 discovery of A. tumefaciens-mediated T-DNA transfer, along with the realization that any DNA of interest can be substituted for oncogenic T-DNA, spawned an entirely new field of plant genetic engineering ¹²¹. In recent years, 484

485 conjugation machines have been repurposed to deliver CRISPR–Cas9 systems to bacterial recipients to cure drug 486 resistance plasmids or kill recipient cells harbouring CRISPR–Cas9 target sequences ¹²²⁻¹²⁵. Very recently, 487 bacterial donors engineered to surface display nanobodies were shown to selectively deliver DNA cargoes to 488 recipient cells displaying the cognate antigens ¹²⁶. These types of translational advances set the stage for 489 deployment of T4SSs for selective killing of bacterial targets or even of cancer cells.

490

491 In the 20 years since the state of knowledge of the fascinatingly versatile T4SSs was first reviewed², the field has 492 made astounding progress in defining many T4SS structures and mechanisms of actions, and identifying the 493 range of cellular consequences accompanying effector translocation. We fully expect the next 20 years to yield 494 even more exciting fundamental and translational advances.

- 495
- 496

497 Acknowledgements

Work in the authors' laboratories was supported by the Welcome Trust grants 215164/Z/18/Z and 217089/Z/19/Z to T.R.D.C. and G.W., respectively, and National Institutes of Health grants NIH 1R35GM131892 and NIH 1R21AI159970 to P.J.C.

501

502 Autor contributions

- T.R.D.C, J.B.P. and P.J.C. wrote the article. T.R.D.C., J.B.P., P.J.C. and K.M. researched data for the article. K.M.
 made figures and table with contributions from T.R.D.C. and J.B.P. All authors reviewed and/or edited the
 manuscript before submission.
- 506

507 Competing interests

508 The authors declare no competing interests.

509

510 Supplementary information

- 511 Supplementary information is available for this paper at https://doi.org/10.1038/s415XX-XXX-XXX-X
- 512

513 References

514

- Alvarez-Martinez, C. E. & Christie, P. J. Biological diversity of prokaryotic type IV secretion systems.
 Microbiology and Molecular Biology Reviews 73, 775-808 (2009).
- 517 2 Cascales, E. & Christie, P. J. The versatile bacterial type IV secretion systems. *Nature Reviews* 518 *Microbiology* **1**, 137-149 (2003). <u>https://doi.org:10.1038/nrmicro753</u>
- 5193Costa, T. R. *et al.* Secretion systems in Gram-negative bacteria: structural and mechanistic insights. *Nat*520*Rev Microbiol* 13, 343-359 (2015). https://doi.org:10.1038/nrmicro3456
- 521 4 Christie, P. J. The Mosaic Type IV Secretion Systems. *EcoSal Plus* **7** (2016). 522 <u>https://doi.org:10.1128/ecosalplus.ESP-0020-2015</u>
- 523 5 Waksman, G. From conjugation to T4S systems in Gram-negative bacteria: a mechanistic biology 524 perspective. *EMBO Rep* **20** (2019). <u>https://doi.org:10.15252/embr.201847012</u>
- 525 6 Cabezon, E., de la Cruz, F. & Arechaga, I. Conjugation Inhibitors and Their Potential Use to Prevent
 526 Dissemination of Antibiotic Resistance Genes in Bacteria. *Front Microbiol* 8, 2329 (2017).
 527 https://doi.org:10.3389/fmicb.2017.02329
- 528 7 Boudaher, E. & Shaffer, C. L. Inhibiting bacterial secretion systems in the fight against antibiotic 529 resistance. *Medchemcomm* **10**, 682-692 (2019). <u>https://doi.org:10.1039/c9md00076c</u>
- 530 8 Christie, P. J., Atmakuri, K., Krishnamoorthy, V., Jakubowski, S. & Cascales, E. Biogenesis, architecture,
 531 and function of bacterial type IV secretion systems. *Annu Rev Microbiol* 59, 451-485 (2005).
 532 https://doi.org:10.1146/annurev.micro.58.030603.123630
- 533 9 Cabezón, E., Ripoll-Rozada, J., Peña, A., de la Cruz, F. & Arechaga, I. Towards an integrated model of 534 bacterial conjugation. *FEMS Microbiol Rev* **39**, 81-95 (2015). <u>https://doi.org:10.1111/1574-6976.12085</u>
- 53510Costa, T. R. D. *et al.* Type IV secretion systems: Advances in structure, function, and activation. *Mol*536*Microbiol* **115**, 436-452 (2021). https://doi.org:10.1111/mmi.14670
- 537 11 Sheedlo, M. J., Ohi, M. D., Lacy, D. B. & Cover, T. L. Molecular architecture of bacterial type IV secretion
 538 systems. *PLoS Pathog* 18, e1010720 (2022). <u>https://doi.org:10.1371/journal.ppat.1010720</u>
- 539
 12
 Sgro, G. G. et al. Bacteria-Killing Type IV Secretion Systems. Front Microbiol 10, 1078 (2019).

 540
 https://doi.org:10.3389/fmicb.2019.01078
- 54113Gonzalez-Rivera, C., Bhatty, M. & Christie, P. J. Mechanism and Function of Type IV Secretion During542Infection of the Human Host. Microbiol Spectr 4 (2016). https://doi.org/10.1128/microbiolspec.VMBF-0024-2015
- 544
 14
 Macé, K. et al. Cryo-EM structure of a type IV secretion system. Nature 607, 191-196 (2022).

 545
 https://doi.org:10.1038/s41586-022-04859-y
- 546The first high-resolution atomic structure of a nearly complete T4SS representing a significant leap547forward in the understanding of type IV secretion, as it reveals crucial details of assembly, function,548and subunit interfaces, opening new possibilities for rational drug design and establishing a workflow549for structural determination of these complex machineries.
- 55015Amin, H., Ilangovan, A. & Costa, T. R. D. Architecture of the outer-membrane core complex from a551conjugative type IV secretion system. Nature Communications 12, 6834 (2021).552https://doi.org:10.1038/s41467-021-27178-8
- 553The first high-resolution structure of the outer membrane core complex from an expanded554conjugative T4SS sheding light on the mechanisms of conjugative pilus outgrowth and DNA555translocation during bacterial conjugation, revealing structural adaptations that contribute to the556dynamic properties of the machinery.
- 55716Chandran, V. *et al.* Structure of the outer membrane complex of a type IV secretion system. Nature558462, 1011-1015 (2009). https://doi.org/10.1038/nature08588

- Fronzes, R. *et al.* Structure of a Type IV Secretion System Core Complex. *Science* 323, 266-268 (2009).
 <u>https://doi.org:doi:10.1126/science.1166101</u>
- 56118Rivera-Calzada, A. *et al.* Structure of a bacterial type IV secretion core complex at subnanometre562resolution.*TheEMBOJournal***32**,1195-1204(2013).563https://doi.org/10.1038/emboj.2013.58
- 564
 19
 Low, H. H. et al. Structure of a type IV secretion system. Nature 508, 550-553 (2014).

 565
 https://doi.org:10.1038/nature13081
- 56620Gordon, J. E. *et al.* Use of chimeric type IV secretion systems to define contributions of outer membrane567subassemblies for contact-dependent translocation. *Mol Microbiol* **105**, 273-293 (2017).568https://doi.org:10.1111/mmi.13700
- 56921Sgro, G. G. *et al.* Cryo-EM structure of the bacteria-killing type IV secretion system core complex from570Xanthomonas citri. Nat Microbiol **3**, 1429-1440 (2018). https://doi.org:10.1038/s41564-018-0262-z
- 571The paper reports the high-resolution cryo-EM structure of a T4SS involved in bacterial killing572advancing our understanding of the structural similarities and differences among functionally distinct573T4SSs.
- 574 22 Cascales, E., Atmakuri, K., Sarkar, M. K. & Christie, P. J. DNA Substrate-Induced Activation of the
 575 Agrobacterium VirB/VirD4 Type IV Secretion System. *Journal of Bacteriology* 195, 2691-2704 (2013).
 576 https://doi.org:doi:10.1128/JB.00114-13
- 577 23 Cascales, E. & Christie, P. J. <i>Agrobacterium</i> VirB10, an ATP energy sensor required for type IV
 578 secretion. *Proceedings of the National Academy of Sciences* 101, 17228-17233 (2004).
 579 https://doi.org:doi:10.1073/pnas.0405843101
- 58024Banta, L. M. *et al.* An Agrobacterium VirB10 Mutation Conferring a Type IV Secretion System Gating581Defect. Journal of Bacteriology 193, 2566-2574 (2011). https://doi.org:doi:10.1128/JB.00038-11
- 58225Darbari, V. C. *et al.* Electrostatic Switching Controls Channel Dynamics of the Sensor Protein VirB10 in583A. tumefaciens Type IV Secretion System. ACS Omega 5, 3271-3281 (2020).584https://doi.org:10.1021/acsomega.9b03313
- Souza, D. P. *et al.* A Component of the Xanthomonadaceae Type IV Secretion System Combines a VirB7
 Motif with a NO Domain Found in Outer Membrane Transport Proteins. *PLOS Pathogens* 7, e1002031
 (2011). <u>https://doi.org:10.1371/journal.ppat.1002031</u>
- Nakano, N., Kubori, T., Kinoshita, M., Imada, K. & Nagai, H. Crystal Structure of Legionella DotD: Insights
 into the Relationship between Type IVB and Type II/III Secretion Systems. *PLOS Pathogens* 6, e1001129
 (2010). <u>https://doi.org:10.1371/journal.ppat.1001129</u>
- 591 28 Lockwood, D. C., Amin, H., Costa, T. R. D. & Schroeder, G. N. The Legionella pneumophila Dot/Icm type 592 IV secretion svstem and its effectors. Microbiology (Readina) 168 (2022). 593 https://doi.org:10.1099/mic.0.001187
- 59429Aly, K. A. & Baron, C. The VirB5 protein localizes to the T-pilus tips in Agrobacterium tumefaciens.595Microbiology 153, 3766-3775 (2007). https://doi.org:10.1099/mic.0.2007/010462-0
- 59630Jakubowski, S. J., Krishnamoorthy, V., Cascales, E. & Christie, P. J. Agrobacterium tumefaciens VirB6597domains direct the ordered export of a DNA substrate through a type IV secretion System. Journal of598molecular biology 341, 961-977 (2004). https://doi.org:10.1016/j.jmb.2004.06.052
- 59931Hospenthal, M. K., Costa, T. R. D. & Waksman, G. A comprehensive guide to pilus biogenesis in Gram-
negative bacteria. Nature Reviews Microbiology 15, 365-379 (2017).601https://doi.org:10.1038/nrmicro.2017.40
- 60232Aly, K. A. & Baron, C. The VirB5 protein localizes to the T-pilus tips in Agrobacterium tumefaciens.603Microbiology 153, 3766-3775 (2007). https://doi.org/10.1099/mic.0.2007/010462-0
- 60433Backert, S., Fronzes, R. & Waksman, G. VirB2 and VirB5 proteins: specialized adhesins in bacterial type-605IV secretion systems? Trends in Microbiology 16, 409-413 (2008).606https://doi.org/10.1016/j.tim.2008.07.001

- 60734Bönig, T., Olbermann, P., Bats, S. H., Fischer, W. & Josenhans, C. Systematic site-directed mutagenesis608of the Helicobacter pylori CagL protein of the Cag type IV secretion system identifies novel functional609domains. Sci Rep 6, 38101 (2016). https://doi.org:10.1038/srep38101
- 610 Pham, K. T. et al. Cagl is an essential component of the Helicobacter pylori Cag type IV secretion system 35 611 PLoS and forms complex with CagL. One 7, e35341 (2012). а 612 https://doi.org:10.1371/journal.pone.0035341
- 613 36 Khara, P., Song, L., Christie, P. J. & Hu, B. In Situ Visualization of the pKM101-Encoded Type IV Secretion
 614 System Reveals a Highly Symmetric ATPase Energy Center. *mBio* 12, e0246521 (2021).
 615 <u>https://doi.org:10.1128/mBio.02465-21</u>
- Liu, X., Khara, P., Baker, M. L., Christie, P. J. & Hu, B. Structure of a type IV secretion system core complex
 encoded by multi-drug resistance F plasmids. *Nature Communications* 13, 379 (2022).
 https://doi.org:10.1038/s41467-022-28058-5
- 61938Kitao, T., Kubori, T. & Nagai, H. Recent advances in structural studies of the Legionella pneumophila620Dot/Icm type IV secretion system. *Microbiol Immunol* 66, 67-74 (2022). https://doi.org/10.1111/1348-0421.12951
- 62239Gomez-Valero, L. *et al.* More than 18,000 effectors in the Legionella genus genome provide multiple,623independent combinations for replication in human cells. *Proc Natl Acad Sci U S A* **116**, 2265-2273624(2019). https://doi.org:10.1073/pnas.1808016116
- 62540Asrat, S., de Jesus, D. A., Hempstead, A. D., Ramabhadran, V. & Isberg, R. R. Bacterial pathogen626manipulation of host membrane trafficking. Annu Rev Cell Dev Biol **30**, 79-109 (2014).627https://doi.org:10.1146/annurev-cellbio-100913-013439
- 62841Durie, C. L. *et al.* Structural analysis of the Legionella pneumophila Dot/Icm type IV secretion system629core complex. *Elife* 9 (2020). https://doi.org:10.7554/eLife.59530
- 630The first high-resolution structure of the Legionella pneumophila Dot/Icm T4SS, which plays a crucial631role in niche establishment and the pathogenesis of Legionnaire's disease.
- 63242Sheedlo, M. J. *et al.* Cryo-EM reveals new species-specific proteins and symmetry elements in the633Legionella pneumophila Dot/Icm T4SS. *eLife* **10**, e70427 (2021). https://doi.org/10.7554/eLife.70427
- 63443Varga, M. G. *et al.* Pathogenic Helicobacter pylori strains translocate DNA and activate TLR9 via the635cancer-associated cag type IV secretion system. Oncogene **35**, 6262-6269 (2016).636https://doi.org:10.1038/onc.2016.158
- 637 44 Tegtmeyer, N., Neddermann, M., Asche, C. I. & Backert, S. Subversion of host kinases: a key network in cellular signaling hijacked by Helicobacter pylori CagA. *Mol Microbiol* 105, 358-372 (2017).
 639 https://doi.org:10.1111/mmi.13707
- 64045Pfannkuch, L. *et al.* ADP heptose, a novel pathogen-associated molecular pattern identified in641Helicobacter pylori. *Faseb j* **33**, 9087-9099 (2019). https://doi.org:10.1096/fj.201802555R
- 642 46 Cover, T. L., Lacy, D. B. & Ohi, M. D. The Helicobacter pylori Cag Type IV Secretion System. *Trends in Microbiology* 28, 682-695 (2020). <u>https://doi.org/10.1016/j.tim.2020.02.004</u>
- 64447Frick-Cheng, A. E. *et al.* Molecular and Structural Analysis of the Helicobacter pylori cag Type IV645Secretion System Core Complex. *mBio* 7, e02001-02015 (2016). https://doi.org/10.1128/mBio.02001-64615
- 64748Sheedlo, M. J. *et al.* Cryo-EM reveals species-specific components within the Helicobacter pylori Cag648type IV secretion system core complex. *Elife* 9 (2020). https://doi.org:10.7554/eLife.59495
- Hu, B., Khara, P. & Christie, P. J. Structural bases for F plasmid conjugation and F pilus biogenesis in
 <i>Escherichia coli</i>>. Proceedings of the National Academy of Sciences 116, 14222-14227 (2019).
 https://doi.org:doi:10.1073/pnas.1904428116
- 652The study presents the first in situ cryo-electron tomography visualization of the entire conjugative F653T4SS, revealing four distinct conformational states of the machinery, and providing a step-by-step

654 mechanism of bacterial conjugation and pilus outgrowth, along with architectural representations of 655 the machinery in each state.

- 656 50 Ghosal, D., Chang, Y. W., Jeong, K. C., Vogel, J. P. & Jensen, G. J. In situ structure of the Legionella
 657 Dot/Icm type IV secretion system by electron cryotomography. *EMBO Rep* 18, 726-732 (2017).
 658 https://doi.org:10.15252/embr.201643598
- 659The first cryo-electron tomography visualization of the Dot/Icm T4SS in Legionella pneumophila660revealing a common overall architecture shared across functionally diverse T4SSs.
- 661 51 Chetrit, D., Hu, B., Christie, P. J., Roy, C. R. & Liu, J. A unique cytoplasmic ATPase complex defines the
 662 Legionella pneumophila type IV secretion channel. *Nature Microbiology* 3, 678-686 (2018).
 663 https://doi.org:10.1038/s41564-018-0165-z
- Park, D., Steiner, S., Shao, M., Roy, C. R. & Liu, J. Developmental Transitions Coordinate Assembly of the
 Coxiella burnetii Dot/Icm Type IV Secretion System. *Infect Immun* 90, e0041022 (2022).
 https://doi.org:10.1128/iai.00410-22
- 53 Newton, H. J., McDonough, J. A. & Roy, C. R. Effector Protein Translocation by the Coxiella burnetii
 bot/Icm Type IV Secretion System Requires Endocytic Maturation of the Pathogen-Occupied Vacuole.
 PLOS ONE 8, e54566 (2013). https://doi.org:10.1371/journal.pone.0054566
- 670 54 Chang, Y. W., Shaffer, C. L., Rettberg, L. A., Ghosal, D. & Jensen, G. J. In Vivo Structures of the
 671 Helicobacter pylori cag Type IV Secretion System. *Cell Rep* 23, 673-681 (2018).
 672 https://doi.org:10.1016/j.celrep.2018.03.085
- 673 55 Hu, B. *et al.* In Situ Molecular Architecture of the Helicobacter pylori Cag Type IV Secretion System.
 674 *mBio* 10 (2019). <u>https://doi.org:10.1128/mBio.00849-19</u>
- 675 56 Chung, J. M. *et al.* Structure of the Helicobacter pylori Cag type IV secretion system. *Elife* 8 (2019).
 676 <u>https://doi.org:10.7554/eLife.47644</u>
- 677 The first high-resolution structure of the *Helicobacter pylori* Cag T4SS, which is required for *H. pylori* 678 infection of the human gastrointestinal tract.
- 57 Tegtmeyer, N. *et al.* Toll-like Receptor 5 Activation by the CagY Repeat Domains of Helicobacter pylori.
 680 *Cell Rep* 32, 108159 (2020). <u>https://doi.org:10.1016/j.celrep.2020.108159</u>
- 68158Audette, G. F., Manchak, J., Beatty, P., Klimke, W. A. & Frost, L. S. Entry exclusion in F-like plasmids682requires intact TraG in the donor that recognizes its cognate TraS in the recipient. *Microbiology*683(*Reading*) 153, 442-451 (2007). https://doi.org:10.1099/mic.0.2006/001917-0
- 684 59 Marrero, J. & Waldor, M. K. Determinants of entry exclusion within Eex and TraG are cytoplasmic. *J* 685 *Bacteriol* **189**, 6469-6473 (2007).
- 686 60 Gillespie, J. J. *et al.* An anomalous type IV secretion system in Rickettsia is evolutionarily conserved.
 687 *PLoS One* 4, e4833 (2009). <u>https://doi.org:10.1371/journal.pone.0004833</u>
- 68861Gillespie, J. J. et al. Phylogenomics reveals a diverse Rickettsiales type IV secretion system. Infect Immun68978, 1809-1823 (2010). https://doi.org/10.1128/iai.01384-09
- 62 Rancès, E., Voronin, D., Tran-Van, V. & Mavingui, P. Genetic and Functional Characterization of the Type
 691 IV Secretion System in <i>Wolbachia</i>. Journal of Bacteriology 190, 5020-5030 (2008).
 692 https://doi.org:doi:10.1128/JB.00377-08
- 69363Nagai, H. & Roy, C. R. The DotA protein from Legionella pneumophila is secreted by a novel process that694requirestheDot/Icmtransporter.Emboj20,5962-5970(2001).695https://doi.org:10.1093/emboj/20.21.5962
- 69664Skoog, E. C. *et al.* CagY-Dependent Regulation of Type IV Secretion in Helicobacter pylori Is Associated697with Alterations in Integrin Binding. *MBio* **9** (2018). https://doi.org:10.1128/mBio.00717-18
- 69865Aras, R. A., Kang, J., Tschumi, A. I., Harasaki, Y. & Blaser, M. J. Extensive repetitive DNA facilitates699prokaryotic genome plasticity. *Proc Natl Acad Sci U S A* **100**, 13579-13584 (2003).700https://doi.org:10.1073/pnas.1735481100

- 70166Barrozo, R. M. *et al.* Functional Plasticity in the Type IV Secretion System of Helicobacter pylori. *PLOS*702*Pathogens* 9, e1003189 (2013). https://doi.org:10.1371/journal.ppat.1003189
- 70367Llosa, M. & Alkorta, I. in Type IV Secretion in Gram-Negative and Gram-Positive Bacteria (eds Steffen704Backert & Elisabeth Grohmann) 143-168 (Springer International Publishing, 2017).
- 70568Gomis-Rüth, F. X. *et al.* The bacterial conjugation protein TrwB resembles ring helicases and F1-ATPase.706Nature **409**, 637-641 (2001). https://doi.org/10.1038/35054586
- 70769Whitaker, N. *et al.* The All-Alpha Domains of Coupling Proteins from the Agrobacterium tumefaciens708VirB/VirD4 and Enterococcus faecalis pCF10-Encoded Type IV Secretion Systems Confer Specificity to709Binding of Cognate DNA Substrates. Journal of Bacteriology 197, 2335-2349 (2015).710https://doi.org:doi:10.1128/JB.00189-15
- 71170Oka, G. U. et al. Structural basis for effector recognition by an antibacterial type IV secretion system.712Proceedings of the National Academy of Sciences 119, e2112529119 (2022).713https://doi.org:doi:10.1073/pnas.2112529119
- Atmakuri, K., Cascales, E. & Christie, P. J. Energetic components VirD4, VirB11 and VirB4 mediate early
 DNA transfer reactions required for bacterial type IV secretion. *Mol Microbiol* 54, 1199-1211 (2004).
 https://doi.org:10.1111/j.1365-2958.2004.04345.x
- 717 72 Ripoll-Rozada, J., Zunzunegui, S., de la Cruz, F., Arechaga, I. & Cabezon, E. Functional interactions of
 718 VirB11 traffic ATPases with VirB4 and VirD4 molecular motors in type IV secretion systems. *J Bacteriol* 719 195, 4195-4201 (2013). https://doi.org:10.1128/JB.00437-13
- 72073Savvides, S. N. *et al.* VirB11 ATPases are dynamic hexameric assemblies: new insights into bacterial type721IV secretion. *EMBO J.* 22, 1969-1980 (2003).
- 72274Hare, S., Bayliss, R., Baron, C. & Waksman, G. A Large Domain Swap in the VirB11 ATPase of Brucella723suis Leaves the Hexameric Assembly Intact. Journal of molecular biology 360, 56-66 (2006).724https://doi.org/10.1016/j.jmb.2006.04.060
- 725 75 Park, D., Chetrit, D., Hu, B., Roy, C. R. & Liu, J. Analysis of Dot/Icm Type IVB Secretion System
 726 Subassemblies by Cryoelectron Tomography Reveals Conformational Changes Induced by DotB Binding.
 727 *mBio* 11, e03328-03319 (2020). <u>https://doi.org:doi:10.1128/mBio.03328-19</u>
- 76 Sagulenko, E., Sagulenko, V., Chen, J. & Christie, P. J. Role of Agrobacterium VirB11 ATPase in T-pilus assembly and substrate selection. *Journal of Bacteriology* 183, 5813-5825 (2001).
 730 https://doi.org:10.1128/JB.183.20.5813-5825.2001
- 731 77 Hilleringmann, M. *et al.* Inhibitors of Helicobacter pylori ATPase Cagα block CagA transport and cag
 732 virulence. *Microbiology* 152, 2919-2930 (2006). <u>https://doi.org/10.1099/mic.0.28984-0</u>
- 73378Nagai, H. et al. A C-terminal translocation signal required for Dot/Icm-dependent delivery of the734Legionella RalF protein to host cells. Proceedings of the National Academy of Sciences 102, 826-831735(2005).
- 73679Cambronne, E. D. & Roy, C. R. The Legionella pneumophila IcmSW Complex Interacts with Multiple737Dot/Icm Effectors to Facilitate Type IV Translocation. PLOS Pathogens 3, e188 (2007).738https://doi.org:10.1371/journal.ppat.0030188
- 73980Kim, H. *et al.* Structural basis for effector protein recognition by the Dot/Icm Type IVB coupling protein740complex. Nature communications 11, 1-11 (2020).
- 741 81 Jeong, K. C., Sutherland, M. C. & Vogel, J. P. Novel export control of a Legionella Dot/Icm substrate is
 742 mediated by dual, independent signal sequences. *Mol Microbiol* 96, 175-188 (2015).
 743 <u>https://doi.org:10.1111/mmi.12928</u>
- 74482Meir, A., Chetrit, D., Liu, L., Roy, C. R. & Waksman, G. Legionella DotM structure reveals a role in effector745recruiting to the Type 4B secretion system. Nature Communications 9, 507 (2018).746https://doi.org:10.1038/s41467-017-02578-x
- 747The paper presents the crystallographic structure of DotM, revealing a novel mechanism of effector748recognition by the Dot/Icm T4SS in Legionella pneumophila.

- 74983Mace, K. et al. Proteins DotY and DotZ modulate the dynamics and localization of the type IVB coupling750complex of Legionella pneumophila. Mol Microbiol 117, 307-319 (2022).751https://doi.org:10.1111/mmi.14847
- 75284Meir, A., Macé, K., Vegunta, Y., Williams, S. M. & Waksman, G. Substrate recruitment mechanism by753gram-negative type III, IV, and VI bacterial injectisomes. Trends Microbiol (2023).754https://doi.org:10.1016/j.tim.2023.03.005
- 75585Kwak, M. J. *et al.* Architecture of the type IV coupling protein complex of Legionella pneumophila. *Nat*756*Microbiol* 2, 17114 (2017). https://doi.org:10.1038/nmicrobiol.2017.114
- 75786Meir, A. *et al.* Mechanism of effector capture and delivery by the type IV secretion system from758Legionella pneumophila. Nature Communications 11, 2864 (2020). https://doi.org:10.1038/s41467-759020-16681-z
- Vincent, C. D., Friedman, J. R., Jeong, K. C., Sutherland, M. C. & Vogel, J. P. Identification of the DotL coupling protein subcomplex of the Legionella Dot/Icm type IV secretion system. *Molecular microbiology* 85, 378-391 (2012). https://doi.org:10.1111/j.1365-2958.2012.08118.x
- 76388Xu, J. et al. Structural insights into the roles of the IcmS–IcmW complex in the type IVb secretion system764of Legionella pneumophila. Proceedings of the National Academy of Sciences 114, 13543-13548 (2017).
- Pattis, I., Weiss, E., Laugks, R., Haas, R. & Fischer, W. The *Helicobacter pylori* CagF protein is a type IV secretion chaperone-like molecule that binds close to the C-terminal secretion signal of the CagA effector protein. *Microbiology* 153, 2896-2909 (2007). https://doi.org:10.1099/mic.0.2007/007385-0
- 76890Wu, X. et al. Mechanism of regulation of the Helicobacter pylori Cagβ ATPase by CagZ. Nature769Communications 14, 479 (2023). https://doi.org:10.1038/s41467-023-36218-4
- de la Cruz, F., Frost, L. S., Meyer, R. J. & Zechner, E. L. Conjugative DNA metabolism in Gram-negative bacteria. *FEMS microbiology reviews* 34, 18-40 (2010). <u>https://doi.org:10.1111/j.1574-6976.2009.00195.x</u>
- 77392Ilangovan, A. *et al.* Cryo-EM Structure of a Relaxase Reveals the Molecular Basis of DNA Unwinding774duringBacterialConjugation.Cell169,708-721.e712(2017).775https://doi.org/10.1016/j.cell.2017.04.010
- 776The study reveals the first structure of the relaxase protein, the core component of the relaxosome777complex involved in DNA processing prior to conjugation, showing that two distinct activities of the778relaxase, the trans-esterase activity required for DNA nicking, and the helicase activity essential for779DNA unwinding, are simultaneously performed by two distinct structural conformers.
- 780 93 Datta, S., Larkin, C. & Schildbach, J. F. Structural insights into single-stranded DNA binding and cleavage
 781 by F factor Tral. *Structure* 11, 1369-1379 (2003).
- 78294Rice, P. A., Yang, S., Mizuuchi, K. & Nash, H. A. Crystal structure of an IHF-DNA complex: a protein-
induced DNA U-turn. *Cell* 87, 1295-1306 (1996).
- 5 Luo, Y., Gao, Q. & Deonier, R. C. Mutational and physical analysis of F plasmid *traY* protein binding to
 oriT. Mol Microbiol 11, 459-469 (1994).
- 78696Lu, J. et al. Structural basis of specific TraD–TraM recognition during F plasmid-mediated bacterial
conjugation. Molecular Microbiology 70, 89-99 (2008). https://doi.org/10.1111/j.1365-2958.2008.06391.x
- 789 97 Beranek, A. et al. Thirty-Eight C-Terminal Amino Acids of the Coupling Protein TraD of the F-Like 790 Conjugative Resistance Plasmid R1 Are Required and Sufficient To Confer Binding to the Substrate 791 Selector Protein TraM. Bacteriology 186, 6999-7006 (2004). Journal of 792 https://doi.org:doi:10.1128/JB.186.20.6999-7006.2004
- 79398Wong, J. J., Lu, J., Edwards, R. A., Frost, L. S. & Glover, J. M. Structural basis of cooperative DNA794recognition by the plasmid conjugation factor, TraM. Nucleic acids research **39**, 6775-6788 (2011).

- Rêgo, A. T., Chandran, V. & Waksman, G. Two-step and one-step secretion mechanisms in Gramnegative bacteria: contrasting the type IV secretion system and the chaperone-usher pathway of pilus
 biogenesis. *Biochem J* 425, 475-488 (2010). https://doi.org:10.1042/bj20091518
- 798100Prevost, M. S. & Waksman, G. X-ray crystal structures of the type IVb secretion system DotB ATPases.799Protein Sci 27, 1464-1475 (2018). https://doi.org/10.1002/pro.3439
- Jakubowski, S. J., Cascales, E., Krishnamoorthy, V. & Christie, P. J. Agrobacterium tumefaciens VirB9, an
 outer-membrane-associated component of a type IV secretion system, regulates substrate selection
 and T-pilus biogenesis. J Bacteriol 187, 3486-3495 (2005).
- 803102Redzej, A. *et al.* Structure of a VirD4 coupling protein bound to a VirB type IV secretion machinery. *The*804*EMBO Journal* **36**, 3080-3095 (2017). https://doi.org/10.15252/embj.201796629
- 805 103 Burns, D. L. Secretion of Pertussis Toxin from Bordetella pertussis. Toxins 13, 574 (2021).
- 806104Dehio, C. & Tsolis, R. M. Type IV Effector Secretion and Subversion of Host Functions by Bartonella and
Brucella Species. *Curr Top Microbiol Immunol* **413**, 269-295 (2017). https://doi.org/10.1007/978-3-319-75241-9 11
- 809105Costa, T. R. D. *et al.* Structure of the Bacterial Sex F Pilus Reveals an Assembly of a Stoichiometric810Protein-PhospholipidComplex.*Cell***166**,1436-1444e1410(2016).811https://doi.org:10.1016/j.cell.2016.08.025
- 812The paper presents the structure of the conjugative F pilus, uncovering the incorporation of813phospholipid molecules within the pilus filament's molecular architecture, establishing the basis for814structural characterization of conjugative pili and leading to a surge in available architectures of other815conjugative pili with their respective phospholipid types.
- 816106Zheng, W. et al. Cryoelectron-Microscopic Structure of the pKpQIL Conjugative Pili from Carbapenem-817ResistantKlebsiellapneumoniae.Structure28,1321-1328.e1322(2020).818https://doi.org:10.1016/j.str.2020.08.010
- 819107Kreida, S. et al. Cryo-EM structure of the Agrobacterium tumefaciens T4SS-associated T-pilus reveals820stoichiometricprotein-phospholipidassembly.Structure(2023).821https://doi.org:10.1016/j.str.2023.02.005
- 822108Amro, J. *et al.* Cryo-EM structure of the Agrobacterium tumefaciens T-pilus reveals the importance of823positive charges in the lumen. *Structure* **31**, 375-384 (2023). https://doi.org/10.1016/j.str.2022.11.007
- Beltran, L. C. *et al.* Archaeal DNA-import apparatus is homologous to bacterial conjugation machinery.
 Nature Communications 14, 666 (2023). <u>https://doi.org:10.1038/s41467-023-36349-8</u>
- Patkowski, J. B. *et al.* The F-pilus biomechanical adaptability accelerates conjugative dissemination of antimicrobial resistance and biofilm formation. *Nature Communications* 14, 1879 (2023).
 https://doi.org:10.1038/s41467-023-37600-y
- Rozwandowicz, M. *et al.* Plasmids carrying antimicrobial resistance genes in Enterobacteriaceae. J
 Antimicrob Chemother **73**, 1121-1137 (2018). https://doi.org:10.1093/jac/dkx488
- 831 112 Bradley, D. E. Morphological and serological relationships of conjugative pili. *Plasmid* 4, 155-169 (1980).
- Paschos, A. *et al.* An in vivo high-throughput screening approach targeting the type IV secretion system component VirB8 identified inhibitors of Brucella abortus 2308 proliferation. *Infection and immunity* 79, 1033-1043 (2011). <u>https://doi.org:10.1128/IAI.00993-10</u>
- 835 114 Shaffer, C. L. *et al.* Peptidomimetic Small Molecules Disrupt Type IV Secretion System Activity in Diverse
 836 Bacterial Pathogens. *MBio* 7, e00221-00216 (2016). <u>https://doi.org:10.1128/mBio.00221-16</u>
- Ripoll-Rozada, J. *et al.* Type IV traffic ATPase TrwD as molecular target to inhibit bacterial conjugation.
 Mol Microbiol 100, 912-921 (2016). <u>https://doi.org:10.1111/mmi.13359</u>
- Casu, B., Arya, T., Bessette, B. & Baron, C. Fragment-based screening identifies novel targets for inhibitors of conjugative transfer of antimicrobial resistance by plasmid pKM101. *Scientific Reports* 7, 14907 (2017). <u>https://doi.org:10.1038/s41598-017-14953-1</u>

- 842 117 Getino, M. & de la Cruz, F. Natural and Artificial Strategies To Control the Conjugative Transmission of
 843 Plasmids. *Microbiol Spectr* 6 (2018). <u>https://doi.org:10.1128/microbiolspec.MTBP-0015-2016</u>
- 844 118 Garcia-Cazorla, Y. *et al.* Conjugation inhibitors compete with palmitic acid for binding to the conjugative traffic ATPase TrwD, providing a mechanism to inhibit bacterial conjugation. *J Biol Chem* 293, 16923-16930 (2018). <u>https://doi.org:10.1074/jbc.RA118.004716</u>
- Arya, T. *et al.* Fragment-based screening identifies inhibitors of ATPase activity and of hexamer formation of Cagalpha from the Helicobacter pylori type IV secretion system. *Sci Rep* 9, 6474 (2019).
 https://doi.org:10.1038/s41598-019-42876-6
- 850120Alvarez-Rodriguez, I. *et al.* Type IV Coupling Proteins as Potential Targets to Control the Dissemination851of Antibiotic Resistance. Front Mol Biosci 7, 201 (2020). https://doi.org/10.3389/fmolb.2020.00201
- Brown, P. J. B., Chang, J. H. & Fuqua, C. Agrobacterium tumefaciens: a Transformative Agent for
 Fundamental Insights into Host-Microbe Interactions, Genome Biology, Chemical Signaling, and Cell
 Biology. J Bacteriol, e0000523 (2023). https://doi.org:10.1128/jb.00005-23
- Hamilton, T. A. *et al.* Efficient inter-species conjugative transfer of a CRISPR nuclease for targeted bacterial killing. *Nature Communications* 10, 4544 (2019). <u>https://doi.org:10.1038/s41467-019-12448-</u>
 3
- Vrancianu, C. O., Popa, L. I., Bleotu, C. & Chifiriuc, M. C. Targeting Plasmids to Limit Acquisition and Transmission of Antimicrobial Resistance. *Front Microbiol* **11**, 761 (2020). https://doi.org:10.3389/fmicb.2020.00761
- Reuter, A. *et al.* Targeted-antibacterial-plasmids (TAPs) combining conjugation and CRISPR/Cas systems
 achieve strain-specific antibacterial activity. *Nucleic Acids Res* 49, 3584-3598 (2021).
 https://doi.org:10.1093/nar/gkab126
- 864 125 Bier, E. & Nizet, V. Driving to Safety: CRISPR-Based Genetic Approaches to Reducing Antibiotic 865 Resistance. *Trends Genet* **37**, 745-757 (2021). <u>https://doi.org:10.1016/j.tig.2021.02.007</u>
- Robledo, M. *et al.* Targeted bacterial conjugation mediated by synthetic cell-to-cell adhesions. *Nucleic Acids Res* 50, 12938-12950 (2022). <u>https://doi.org:10.1093/nar/gkac1164</u>
- 868
- 869

870 Figure 1. The functional versatility of type IV secretion systems.

871 Various pathogenic bacteria and symbionts deploy type IV secretion systems (T4SSs) to deliver effector proteins, 872 DNA-protein complexes or other macromolecules into eukaryotic or protozoan host cells. a The T4SSs 873 establishes contact dependent interkingdom interactions by injecting effectors directly into eukaryotic cells to 874 promote bacterial intracellular survival and symbiosis. b | Many bacterial species and a few Archaea deploy a 875 contact-dependent T4SSs for the delivery of DNA and toxins to other bacteria or Archaea. Various species in the 876 Xanthomonadales instead deploy T4SSs for the contact-dependent delivery of protein toxins to kill other 877 bacteria for niche establishment. c | Some bacteria can deploy T4SSs for the contact-independent uptake or 878 release of DNA.

879

880 Figure 2. Structure of minimal type IV secretion system and pilus biogenesis mechanism model.

881 a | Global organisation to atomic details of the minimised R388-encoded type IV secretion system (T4SS). The

- 882 entire cryo-electron microscopy (cryo-EM) structure of the R388 T4SS is shown with half-left in cryo-EM density
- colored by sub-complexes (Electron Microscopy Data Bank (EMDB) entry 12707, EMDB 12708, EMDB 12709,
- 884 EMDB 13767 and EMDB 12933) and half-right in ribbon and surface semi-transparent representation colored by

885 proteins (Protein Data Bank (PDB) identifier 703J, PDB ID 703T, PDB ID 703V, PDB ID 7Q1V and PDB ID 70IU). 886 At top-left corner, the cryo-ET of pKM101 T4SS density (EMDB 24098 and 24100) colored by sub-complexes is 887 displayed to show that the structure of purified R388 T4SS is similar to the in situ T4SS structure. For each sub-888 complex, structure details, symmetry and membrane localisation are indicated. Black dashed lines demarcate 889 the boundaries of the outer membrane and inner membrane. **b** | Model of pilus biogenesis mechanism. The 890 T4SS is schematically represented in slice-view and colored by protein. Four states are shown: T4SS in similar 891 state as observed by cryo-EM and shown in part a (1). The pilus biogenesis state with VirB11 bound at the 892 bottom of VirB4; VirB2 is extracted from the inner membrane and recruited to VirB6 through the coordinated 893 actions of the VirB4–VirB11 ATPases (2). As layers of VirB2 are recruited, the pilus grows from the VirB6 894 assembly sites and VirB5 remains at the pilus tip (3). As the pilus grows, the O-layer of the OMCC opens up and 895 the pilus with VirB5 at the tip extends into the extracellular milieu to establish contact with potential recipient 896 cells (4). Part a is adapted with permission from Ref 14.

897

898 Figure 3. Structural organisation of expanded type IV secretion system.

899 a | Three expanded type IV secretion system (T4SS) outer membrane core complex (OMCC) structures are 900 shown. The OMCC of the F-plasmid (Protein Data Bank (PDB) identifier 70KN and PDB ID 70KO), Legionella 901 pneumophila (PDB ID 7MUS) and Helicobacter pylori (PDB ID 6X6S and 6X6J) are shown in surface representation 902 and colored in dark red, blue and green for the VirB7-, VirB9- and VirB10-like proteins, respectively, and in grey 903 for other components. Notably, L. pneumophila and H. pylori OMCCs contain an outer membrane cap (OMC) 904 and a periplasmic ring (PR). b | Cryo-tomography maps of the F-plasmid with and without pilus. The maps 905 (Electron Microscopy Data Bank (EMDB) entry 9344 and 9347) are colored by subcomplexes (that is; green for 906 the OMCC, red for the Stalk, yellow for the Arches, blue for the inner membrane complex (IMC) and grey for the 907 pilus. The junction of the pilus and Stalk is not well defined. c | Cryo-electron microscopy (cryo-ET) map of the 908 L. pneumophila T4SS (EMDB entry 7611 and 7612) colored as in b, in front and slice view. d | Cryo-ET of the H. 909 pylori T4SS (EMDB entry 0634 and 0635) colored as in b, in front and slice view.

910

911 Figure 4. Examples of type IV secretion system subunit adaptations for functional diversification.

912 The R388-encoded type IV secretion systemT4SS is shown at left for reference. VirB5 subunits are deployed for 913 binding of target-cell receptors; these subunits can localize at the tips of conjugative pili or on the bacterial cell 914 surface. Some bacteria encode several copies of VirB2 or VirB5 subunits whose variable sequences are 915 postulated to bind different target-cell receptors or contribute to evasion of the host immune system. Extended 916 VirB6 carry large hydrophilic domains, several of which have been shown or are implicated in localizing at the 917 cell surface to promote adhesion or immunomodulation, or blocking redundant plasmid transfer. F systems 918 elaborate F pili that dynamically extend and retract to establish contacts with potential recipient cells at a 919 distance. F systems also code for TraN subunits, whose extracellular domains interact with outer membrane 920 proteins (OMPs) on recipient cells to promote F plasmid transfer and specify plasmid host range. Several T4SSs 921 possess variant forms of the VirB7–VirB9–VirB10 core complex subunits, as exemplified for CagY in the H. pylori 922 Cag system. The H. pylori system elaborates a conjugative pilus, which is decorated by other Cag subunits as well

923 as the CagA secretion substrate. Various T4SSs functionally interact with other surface adhesins, such as
 924 pKM101-encoded Pep or *H. pylori* OMPs, to promote target-cell binding.

925

926 Figure 5. Structure of VirD4- and VirB11-like ATPases.

927 a | Side and top view of two VirD4-like ATPase structures: TrwB from R388 plasmid (Protein Data Bank (PDB) 928 identifier (ID) 1GKI) and Cagß from *H. pylori* (PDB ID 8DOL). Structures are shown in ribbon representation and 929 colored by monomer. **b** | Side and top view of two VirB11-like ATPase structures; DotB from *L. pneumophila* 930 (PDB ID 6GEF) and Cagα from *H. pylori* (PDB ID 1NLZ). Structures are shown in ribbon representation and colored 931 by monomer. $\mathbf{c} \mid$ Organisation and structure of type 4 coupling complex (T4CC) from L. pneumophila. The 932 monomeric cryo- electron microscopy (cryo-EM) structure of the T4CC (PDB ID 6SZ9) is shown in surface 933 representation and colored by protein: VirD4-like DotL in red), DotM in cyan, DotN in blue, DotZ in yellow, DotY 934 in green, IcmS in pink and IcmW in purple. The module made up of the extreme C-terminus of (DotL_{C-ter}) and 935 IcmSW is flexible, and its motion area is represented in mesh. Crystal structure of the DotL_{C-ter}-IcmSW module 936 in presence of LvgA adaptor (orange) and VpdB substrate (black) is shown (PDB ID 7BWK). On the right, a model 937 of the hexameric T4CC structure is shown in front and top views using the same color coding as in part c.

938

939 Figure 6. Models for substrates recruitment and transport through the type IV secretion system.

940 a | Conjugative type IV secretion system (T4SS) recruitment and secretion mechanism model. Silhouette of F-941 plasmid T4SS is shown in green. First, the DNA is processed by a relaxosomal complex made of TraM (green), 942 TraY (turquoise), IHF (purple) and TraI (light brown). The relaxosome is recruited by TraD_{VirD4} ATPase, which 943 energizes the secretion of the Tral-ssDNA through the T4SS apparatus into the host. The relative position of 944 TraD_{VirD4} and the global organisation of the inner membrane complex (IMC) during DNA secretion is unknown. **b** 945 | T4SS effector recruitment and secretion mechanism model. Silhouette of L. pneumophila Dot/Icm T4SS is 946 shown in blue. The type IV coupling complex (T4CC) acts as a effector recruitment platform and is schematically 947 represented and positioned beside the complex formed by the hexameric dimers of DotO and DotB, although 948 its precise localisation is unknown. Effector proteins are captured by the T4CC at different binding sites and 949 DotL_{VirD4} energizes substrate translocation via one of two possible routes across the inner membrane. Route 1: 950 The T4CC feeds substrates into the DotO–DotB energy centre at the base of the T4SS channel for transit in one-951 step across the entire cell envelope. Route 2: the T4CC feeds substrates into the lumen of the DotL hexamer for 952 delivery across the inner membrane. In a second translocation step, substrates are recruited from the periplasm 953 by the T4SS channel for passage to the cell surface and into target cells. Part b is adapted with permission from 954 Ref. 28.

955

956 Figure 7. Structure comparison between minimal, expended and archaea pilus.

957 a | Side view of all known pilus structures. Pilus structures (Protein Data Bank (PDB) identifier (ID) 8EXH
 958 Agrobacterium tumefaciens T-pilus, PDB ID 8CW4 Escherichia coli N-pilus; PDB ID 5LER E. coli F-pilus; PDB ID
 959 5LEG Salmonella typhimurium pED208; PDB ID 7JSV Klebsiella pneumoniae pKpQIL; PDB ID 8DFU Aeropyrum
 960 pernix CedA1 and PDB ID 8DFT Pyrobaculum calidifontis TedC) are in surface representation with one strand

- 961 colored in grey. **b** | Top view of pilus structures. Diameter and lumen sizes are indicated. **c** | For each pilus, one
- 962 monomer of VirB2 with its lipid is shown in ribbon representation. Arrows in the minimised type IV secretion
- 963 systems (T4SSs) highlight the presence of a 'kink', which is characteristic of this group. Parts a and b adapted
- 964 with permission from Ref. 109.

Fig. 1: The functional versatility of type IV secretion systems.

From: Structural and functional diversity of type IV secretion systems

a

Contact-dependent interkingdom



b

Contact-dependent interbacterial

Conjugation

Gram-positive and Gram-negative bacteria

Xanthomonadales

С

Contact-independent DNA uptake or release

DNA uptake

Helicobacter pylori

DNA release and biofilm formation

Neisseria gonorrhoeae

Various pathogenic bacteria and symbionts deploy type IV secretion systems (T4SSs) to deliver effector proteins, DNA-protein complexes or other macromolecules into eukaryotic or protozoan host cells. **a**, The T4SS establishes contact-dependent interkingdom interactions by injecting effectors directly into eukaryotic cells to promote bacterial intracellular survival and symbiosis. **b**, Many bacterial species and a few Archaea deploy a contactdependent T4SS for the delivery of DNA and toxins to other bacteria or Archaea. Various species in the *Xanthomonadales* instead deploy T4SSs for the contact-dependent delivery of protein toxins to kill other bacteria for niche establishment. **c**, Some bacteria can deploy T4SSs for the contact-independent uptake or release of DNA. ssDNA, single-stranded DNA.

Fig. 2: Structure of minimal type IV secretion system and pilus biogenesis mechanism model.

From: Structural and functional diversity of type IV secretion systems

a, Global organization to atomic details of the minimized R388-encoded type IV secretion system (T4SS). The entire cryo-electron microscopy (cryo-EM) structure of the R388 T4SS is shown with half-left in cryo-EM density coloured by sub-complexes (Electron Microscopy Data Bank (EMDB) entry 12707, EMDB 12708, EMDB 12709, EMDB 13767 and EMDB 12933) and half-right in ribbon and surface semi-transparent representation coloured by proteins (Protein Data Bank (PDB) identifier (ID) 7O3J, PDB ID 7O3T, PDB ID 7O3V, PDB ID 7Q1V and PDB ID 7O1U). In the top-left corner, the cryo-electron tomography (cryo-ET) of pKM101 T4SS density (EMDB 24098 and 24100) coloured by sub-complexes is displayed to show that the structure of purified R388 T4SS is similar to the in situ T4SS structure. For each sub-complex, structure details, symmetry and membrane localization are indicated. Black dashed lines demarcate the boundaries of the outer membrane and inner membrane. **b**, Model of pilus biogenesis mechanism. The T4SS is schematically represented in slice view and coloured by protein. Four states are shown: (1) T4SS in similar state to that observed by cryo-EM and shown in part **a**. (2) The pilus biogenesis state with VirB11 bound at the bottom of VirB4; VirB2 is extracted from the inner membrane and recruited to VirB6 through the coordinated actions of the VirB4–VirB11 ATPases. (3) As layers of VirB2 are recruited, the pilus grows from the VirB6 assembly sites and VirB5 remains at the pilus tip. (4) As the pilus grows, the O-layer of the outer membrane core complex (OMCC) opens up and the pilus with VirB5 at the tip extends into the extracellular milieu to establish contact with potential recipient cells. IMC, inner membrane complex. Part **a** adapted from ref. ¹⁴, CC BY 4.0 (https://creativecommons.org/licenses/by/4.0/).

Fig. 3: Structural organization of expanded type IV secretion system.

From: Structural and functional diversity of type IV secretion systems

C L. pneumophila cryo-ET

d H. pylori cryo-ET

a, Three expanded type IV secretion system (T4SS) outer membrane core complex (OMCC) structures are shown. The OMCCs of the F plasmid (Protein Data Bank (PDB) identifier (ID) 7OKN and PDB ID 7OKO), *Legionella pneumophila* (PDB ID 7MUS) and *Helicobacter pylori* (PDB ID 6X6S and 6X6J) are shown in surface representation and coloured in dark red, blue and green for the VirB7-like, VirB9-like and VirB10-like proteins, respectively, and in grey for other components. Notably, *L. pneumophila* and *H. pylori* OMCCs contain an outer membrane cap (OMC) and a periplasmic ring (PR). **b**, Cryo-electron tomography (cryo-ET) maps of the F plasmid with and without pilus. The maps (Electron Microscopy Data Bank (EMDB) entry 9344 and 9347) are coloured

Fig. 4: Examples of type IV secretion system subunit adaptations for functional diversification.

From: Structural and functional diversity of type IV secretion systems

The R388-encoded type IV secretion system (T4SS) is shown on the left for reference. VirB5 subunits are deployed for binding of target-cell receptors; these subunits can localize at the tips of conjugative pili or on the bacterial cell surface. Some bacteria encode several copies of VirB2 or VirB5 subunits whose variable sequences are postulated to bind different target-cell receptors or contribute to evasion of the host immune system. Extended VirB6 subunits carry large hydrophilic domains, several of which have been shown or are implicated in localizing at the cell surface to promote adhesion or immunomodulation, or blocking redundant plasmid transfer. F systems elaborate F pili that dynamically extend and retract to establish contacts with potential recipient cells at a distance. F systems also code for TraN subunits, whose extracellular domains interact with outer membrane proteins (OMPs) on recipient cells to promote F plasmid transfer and specify plasmid host range. Several T4SSs possess variant forms of the VirB7-VirB9-VirB10 core complex subunits, as exemplified for CagY in the Helicobacter pylori Cag system. The H. pylori system elaborates a conjugative pilus, which is decorated by other Cag subunits and the CagA secretion substrate. Various T4SSs functionally interact with other surface adhesins, such as pKM101-encoded Pep or H. pylori OMPs, to promote target-cell binding.

Expanded T4SS

Fig. 5: Structure of VirD4-like and VirB11-like ATPases.

From: Structural and functional diversity of type IV secretion systems

a, Side view (left) and top view (right) of two VirD4-like ATPase structures: hexameric TrwB from R388 plasmid (Protein Data Bank (PDB) identifier (ID) 1GKI) and hexameric Cagβ from *Helicobacter pylori* (PDB ID 8DOL). Structures are shown in ribbon representation and coloured by monomer. **b**, Side view (left) and top view (right) of two VirB11-like ATPase structures; hexameric DotB from *Legionella pneumophila* (PDB ID 6GEF) and hexameric Cagα from *H. pylori* (PDB ID 1NLZ). Structures are shown in ribbon representation and coloured by monomer. **c**, Organization and structure of type IV coupling complex (T4CC) from *L. pneumophila*. The monomeric cryo-electron microscopy (cryo-EM) structure of the T4CC (PDB ID 6SZ9) is shown in surface representation and coloured by protein: VirD4-like DotL in red, DotM in cyan, DotN in blue, DotZ in yellow, DotY in green, IcmS in pink and IcmW in purple. The module made up of the extreme C terminus of (DotL_{C-ter}) and IcmSW is flexible, and its motion area is represented in mesh. Crystal structure of the DotL_{C-ter}-IcmSW module in the presence of LvgA adaptor (orange) and VpdB substrate (black) is shown (PDB ID 7BWK). On the right, a model of the hexameric T4CC structure is shown in top view using the same colour coding as in part **c**.

Fig. 6: Models for substrate recruitment and transport through the type IV secretion system.

From: Structural and functional diversity of type IV secretion systems

a, Conjugative type IV secretion system (T4SS) recruitment and secretion mechanism model. Silhouette of F plasmid T4SS is shown in blue. First, the DNA is processed by a relaxosomal complex made of TraM (dark green), TraY (turquoise), IHF (purple and dark blue) and TraI (light green). The relaxosome is recruited by TraD_{VirD4} ATPase, which energizes the secretion of the TraI–single-stranded DNA (ssDNA) through the T4SS apparatus into the host. The relative position of TraD_{VirD4} and the global organization of the inner membrane complex during DNA secretion are unknown. **b**, T4SS effector recruitment and secretion mechanism model. Silhouette of *Legionella pneumophila* Dot/Icm T4SS is shown in blue. The type IV coupling complex (T4CC) acts as an effector recruitment platform and is schematically represented and positioned beside the complex formed by the hexameric dimers of DotO and DotB,

Fig. 7: Structure comparison between minimal, expanded and archaea pilus.

From: Structural and functional diversity of type IV secretion systems

a, Side view of all known pilus structures. Pilus structures (Protein Data Bank (PDB) identifier (ID) 8EXH *Agrobacterium tumefaciens* T pilus, PDB ID 8CW4 *Escherichia coli* N pilus, PDB ID 5LER *E. coli* F pilus, PDB ID 5LEG *Salmonella enterica* subsp. *enterica* serovar Typhimurium pED208, PDB ID 7JSV *Klebsiella pneumoniae* pKpQIL, PDB ID 8DFU *Aeropyrum pernix* CedA1 and PDB ID 8DFT *Pyrobaculum calidifontis* TedC) are in surface representation with one strand coloured in grey. **b**, Top view of pilus structures. Diameter and lumen sizes are indicated. **c**, For each pilus, one monomer of VirB2 with its lipid is shown in ribbon representation. The arrows in the minimized type IV secretion systems (T4SSs) highlight the presence of a 'kink', which is characteristic of this group. Parts **a** and **b** adapted from ref. ¹⁰⁸, CC BY 4.0 (https://creativecommons.org/licenses/by/4.0/).

966	Table of content:
967	In this Review, Costa and colleagues summarise the current knowledge of type IV secretion system (T4SS)
968	functioning in Gram-negative bacteria, with a focus on their architectures and adaptations for specialised
969	functions. They also explore the biogenesis pathways and spatial localization of T4SSs.
970	
971	Glossary
972	
973	Conjugative Pili:
974	Helical hair-like appendages formed by protein-phospholipid complexes that assemble on the surface of bacteria
975	and can act as conduits for DNA transfer between donor and recipient bacteria.
976	
977	Biotic and abiotic adhesion:
978	The attachment of microorganisms to living (biotic) or non-living (abiotic) surfaces, a process that typically
979	facilitates biofilm formation, niche establishment, or infection.
980	
981	Bacterial conjugation:
982	Type of horizontal gene transfer in bacteria where genetic material, such as plasmids containing genes for
983	antibiotic resistance, is transferred from a donor bacterium to a recipient bacterium.
984	
985	Effector:
986	Bacterial protein, often secreted through a dedicated secretion system, that interacts with and manipulates
987	cellular processes within a host organism, promoting bacterial survival, colonization, or infection.
988	
989	Cryo-electron microscopy (cryo–EM):
990	Electron microscopy imaging technique that involves freezing samples in vitreous ice to preserve their native
991	state and is used to visualize the three-dimensional structure of biological molecules and complexes at near-
992	atomic resolution.
993	
994	Conjugative plasmid:
995	A type of a bacterial plasmid that encodes a conjugative machinery, through which the plasmid and its cargoes
996	of antimicrobial resistance genes, virulence factors, or other fitness traits are delivered between bacterial cells.
997	
998	Polytopic topology:

999	Protein structure that contains multiple transmembrane segments embedded in the cell membrane.
1000	
1001	Co-evolution of proteins:
1002	The reciprocal influence and evolution patterns between two proteins that interact or are dependent on each
1003	other for function, and have evolved in a coordinated, non-random manner reflecting their mutual adaptation
1004	over time.
1005	
1006	Site-directed mutagenesis:
1007	Molecular biology technique used to introduce specific nucleotide mutations into DNA sequences, with the
1008	purpose of studying their effects on protein structure and function.
1009	
1010	Cryo-electron tomography (cryo–ET):
1011	Specialized variation of cryo-EM that enables the visualization of large cellular components or organelles within
1012	their cellular environment.
1013	
1014	Cryo-focused ion beam (cryo-FIB):
1015	A technique used to prepare samples for cryo-electron microscopy by thinning frozen samples with a focused
1016	ion beam, leading to an improved signal-to-noise ratio and resolution in the imaging of biological samples.
1017	
1018	Toll-like receptor:
1019	A family of pattern recognition receptors (PRR) in the immune system that specifically recognize conserved
1020	patterns in pathogens and triggers an immune response.
1021	
1022	Crosslinking:
1023	Artificial formation of covalent bonds by a crosslinker between different molecules that interact or co-localize
1024	within a biological sample, with the common application of studying protein-protein or protein-ligand
1025	interactions.
1026	
1027	Rotational raise:
1028	Angle at which adjacent rings stack in helical assemblies, influencing the overall helical symmetry and packing
1029	of the structure.
1030	
1031	Translocation signal:
1032	Specific amino acid sequence that confers recognition of a protein as a substrate for a dedicated transport
1033	machinery for delivery to a specific cellular location, the extracellular milieu, or another bacterial or eukaryotic
1034	cell.
1035	
1036	Correlative light and electron microscopy (CLEM):

1037	Imaging technique that combines fluorescence microscopy and electron microscopy to correlate high-resolution
1038	structural information with specific molecular or cellular labelling in the same sample.
1039	
1040	Nanobodies:
1041	Single-domain antibody fragments derived from heavy-chain-only IgG antibodies that are naturally found in the
1042	Camelidae family, which includes camels, llamas, and alpacas.
1043	
1044	Adhesins:
1045	Proteins found on the surface of cells that facilitate attachment to other biotic or abiotic surfaces.
1046	
1047	Biofilm:
1048	An assemblage of bacteria on a biotic or abiotic surface, often with a defined architecture, that is embedded in
1049	an extracellular matrix typically composed of proteins, DNA, lipids and other biological molecules.
1050	
1051	Relaxosome:
1052	A complex of proteins responsible for specific nicking of the double-stranded DNA, unwinding of DNA strands
1053	and delivering the single-stranded DNA transfer intermediate to the type IV secretion apparatus prior to
1054	conjugation.
1055	