

### A digital real-time high-resolution imaging system to quantify and identify all emitted charged particles

Arthur Bongrand, Samuel Duval, Jérôme Donnard, Julie Champion

#### ▶ To cite this version:

Arthur Bongrand, Samuel Duval, Jérôme Donnard, Julie Champion. A digital real-time high-resolution imaging system to quantify and identify all emitted charged particles. 2023 IEEE Nuclear Science Symposium and Medical Imaging Conference, Nov 2023, Vancouver, Canada. hal-04404033

HAL Id: hal-04404033

https://hal.science/hal-04404033

Submitted on 18 Jan 2024

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

# A digital real-time high-resolution imaging system to quantify and identify all emitted charged particles








A. BONGRAND<sup>1</sup>, S. DUVAL<sup>2</sup>, J. DONNARD<sup>2</sup>, J. CHAMPION<sup>1</sup>

1: Laboratoire Subatech, IMT Atlantique, Nantes Université, CNRS/IN2P3, France












2: S.A.S. Ai4R, France arthur.bongrand@subatech.in2p3.fr











Currently, the characterization of the nature and spatial distribution of radionuclides in a sample is time consuming and fastidious. Indeed, it necessarily requires the use of two distinct analytical techniques and detectors of different nature.

Moreover, the spatial resolutions of the two methods are often detector dependent, which makes it difficult to integrate and interpret these two measurements together.

Can we perform a fast and accurate quality control of radiopharmaceuticals with a single system?






### **Digital Autoradiography...**


Allows to study and locate in 2 dimensions the radioactivity emanating from a surface sample at scales until 10 micrometer.

→ an accurate up-to-date emission imaging technique!


- Well suited for pre-clinical TAT: allows understanding and quantifying biodistributions in tissues at high spatial resolution and low activity.
- **Digital** autoradiography → Ability to analyze and reanalyze endlessly the digital data acquired!
- Variety of applications: Oncology, Neuroscience, Radio-TLC, Geological Studies, Radio Toxicology...



### BeaQuant: a particle by particle digital autoradiography system




### The Spatio Temporal Spectroscopic Autoradiography Analysis (STSAA)



→ The concentration of daughter nuclei in the radioactive equilibrium can be modeled at any time using the **Bateman equation** provided that the characteristics of the decay chain and if the initial concentrations are known.

→ For our application, as these initial concentrations are unknown, we have developed a method to calculate these concentrations from the evolution of the measured count rate.




# Results

α – 4.2 μs


<sup>209</sup>Pb

<sup>209</sup>Bi

<sup>225</sup>Ac PSMA 3.62kBq



→ Usually the quality controls are made from a 1D projection where we can distinguish the free <sup>225</sup>Ac which remains on the deposit spot from the labelled <sup>225</sup>Ac which migrates. To limit interpretation errors, the analysis must be performed at **secular equilibrium** (~7h).



of the decay chain X Secular Out of secular equilibrium (1h)

 $A_1(1h) = 17,13$ 


Ratio = 4.91%


 $A_2(1h) = 348,28$  VS


Digital autoradiography:

Very promising! ©

→ The STSAA reduces the time needed for a quality control after elution to only one hour allowing both isotope localization and identification. Thus, the purity is deduced from the ratio of free to labeled <sup>225</sup>Ac.







 $A_1(7h) = 16,78$ 

 $A_2(7h) = 346,76$ Ratio = 4,83%


### **Discussion**

**M** 

Digital autoradiography combined with a temporal spectroscopy method allows the quick and accurate identification of different radioisotopes and to observe their localization.

The increased robustness of this method comes from its low dependence on the operator's choices compared to 1D method.

From the results, we can observe the **contribution** of a radioisotope called "**free** <sup>213</sup>**Bi**". Clearly visible on the migration front, free <sup>213</sup>Bi corresponds to <sup>213</sup>Bi which is not radiolabeled and whose **parent**, <sup>217</sup>At, has **not decayed locally**. The most probable hypothesis is that the **free** <sup>213</sup>**Bi** comes from the mother solution. This free <sup>213</sup>Bi has also been observed in biological samples. Studies are underway to understand its biological impact.



### Conclusions


## Digital Autoradiography with spectroscopic method:

**Promising** in **acquisition** of **nuclear** data for the production  $\rightarrow$  Develop standard methods for quality control. for nuclear medicine applications like quality control radionuclides of radiopharmaceuticals!

A precise tool to understand the temporal behavior of a chain of labeled

More experiments!

Perspectives







Precise image & quantification from mBq to kBq.

and Associated Equipment, vol. 1035, p. 166807, 2022.

### References

isotopes.

[1] N. Gillings et al., EANM guideline on the validation of analytical methods for radiopharmaceuticals, EJNMMI Radiopharmacy and Chemistry, 2020.

[2] C.B.L. Victorio et al., TSPO expression in a Zika virus murine infection model as an imaging target for acute

infection-induced neuroinflammation, Eur J Nucl Med Mol Imaging, p. 742–755, 2023. [3] H. Lefeuvre et al., « Spectroscopic autoradiography of alpha particles using a parallel ionization multiplier gaseous detector », Nuclear Instruments and Methods in Physics Research Section A: Accelerators, Spectrometers, Detectors

### Acknowledgments

We thank the GIP Arronax for providing of their facilities and a <sup>225</sup>Ac-PSMA source, especially Dr. Cyrille Alliot for his excellent expertise in the field of radiopharmaceutical quality control.



From science to health



2023 IEEE Nuclear Science Symposium (NSS), Medical Imaging Conference (MIC) & Room Temperature Semiconductor Detector (RTSD)